

RealSense Physiotherapy Guidance System (RPGS)

A home-based physiotherapy aid

Chua Wei Pin, He Cheng Hui, Nabilah Binti Abdul Raub
Engineering Science Programme, National University of Singapore

1. Abstract

The RealSense Physiotherapy Guidance System (RPGS) is a home-based physiotherapy aid that provides real time sensing and guidance for physiotherapy patients to complete their exercises correctly and most importantly, safely. It consists of three synergistic components: an Intel RealSense D435 3D camera, an wearable module consisting of an Arduino UNO and an MPU-6050 Six-Axis (Gyro + Accelerometer) MEMS MotionTracking™ Device, and an interactive GUI that provides feedback to the patient in real time. This report will discuss how the RPGS was engineered to provide the best experience and accurate posture detection. It will then discuss the usage of RPGS in supporting a proof-of-concept exercise, the Single Arm Raise and Rotation. The RPGS was able to detect errors in posture and identify specific errors for this exercise.

2. Keywords

Physiotherapist, exercise, posture, RealSense, OpenPose, Arduino, Accelerometer, deprojection

3. Current Problem

Physiotherapy¹ helps restore movement and function when someone is affected by injury, illness or disability. It is an important aspect of our lives because common injuries like a strained shoulder or back can easily be treated with physiotherapy. The Global Burden of Disease Study conducted in 2015 states that 74% of years living with disability in the world is closely related to health conditions where rehabilitation is beneficial. Physiotherapy is one of these possible rehabilitation methods. In physiotherapy treatment, patients see a physiotherapist for a checkup, and to learn and practice rehabilitation exercises from the physiotherapist. The patient is then supposed to head home to continue practicing those exercises for continued and consistent rehabilitation, until the next physiotherapy session. This is important² because maintaining the continued doing of exercises facilitates the rate of recovery and helps to build physical strength. However, without the supervision of a physiotherapist when performing exercises

by themselves at home, patients have a higher chance of performing the exercises incorrectly³. This will very likely slow down progress or even be detrimental to the recovery process. In more extreme cases, it may even lead to injuries^{4,5}; all of which are undesirable effects going against the very reason for physiotherapy itself.

Evidently, these undesirable effects may only manifest when a physiotherapist is not present to provide guidance to the patient and to correct their posture/actions.

Hence, we have decided to design the **RealSense Physiotherapy Guidance System (RPGS)**, a real-time sensing and interactive GUI designed to provide guidance to patients performing physiotherapy exercises at home. The RPGS provides guidance by warning the patient when they are performing the exercises incorrectly, and subsequently teaches them how to correct their incorrect posture/pose.

4. Product Overview

The RPGS is made up of three synergistic components: an Intel RealSense D435 3D camera⁶, a wearable module consisting of an Arduino UNO⁷ and an MPU-6050 Six-Axis (Gyro + Accelerometer) MEMS MotionTracking™ Device⁸, and an interactive GUI that provides feedback to the patient in real time.

4.1 Intel RealSense D435

The Intel RealSense D435 (Figure 1) is a 3D camera that provides quality depth imaging of up to 10m through a pair of depth sensors that leverage on stereo imaging. It has a wide field of view in both the vertical and horizontal dimensions, and a global shutter that allows for accurate motion capture. The D435 will be connected to a laptop to allow the software component of the RPGS to process the live video feed. Further elaboration on the D435's role in RPGS will be in Section 7 along with discussing the core architecture.



Figure 1: Intel RealSense D435 camera with tripod

4.2 Wearable Module

The wearable module consists of an MPU6050 Accelerometer connected to an Arduino UNO. The entire set-up is contained within a 3D printed box (Figure 2) with velcro straps attached at the bottom. The velcro straps can be tightened around the user's limbs to detect for rotation. Figure 3 shows the patient wearing the module around their wrist, with the correct orientation being the USB cable comparatively nearer to the hand. A hole in the casing allows for a wired connection via a 4m-long USB-B cable to the laptop for data processing.




Figure 2: Contents within the casing



Figure 3: Wearable module strapped onto patient's wrist

According to how the accelerometer was aligned within the casing, the coordinate system of the entire wearable module is shown in Figure 4. Acceleration data will be recorded and normalized values for all 3 axes can allow for the determination of the direction of acceleration.




Figure 4: Cartesian coordinate system of the wearable module

4.3 Interactive GUI

An interactive GUI was created using the PyQt⁹ platform and designed to be user-friendly, customisable to the patient's needs and abilities, and responsive in terms of feedback. The interactive GUI consists of two modes - the patient mode and the physiotherapist mode. The patient mode informs the patient on how to make corrections to the exercises whereas the physiotherapist mode allows for customisation of exercises for the patient by adjusting its level of difficulty.

Figure 5 shows the login page for the RPGS. Patients will have their own personal RPGS account with data concerning their exercises stored. Privacy is ensured through this personalised login function. The physiotherapist will have their own account as well, enabling them to take on an administrator role. This allows the physiotherapist to be able to monitor the progress of the patient and customise exercises to the patient's needs and abilities.




Figure 5: login page of the RPGS

4.3.1 Patient Mode

Upon logging into RPGS in patient mode, the patient is brought to the following screen shown in Figure 6. They are able to select an exercise from a list on the left and view the video showing the exercise, on the right. The video consists of a step-by-step procedure on how to perform the exercise accurately to prevent errors in posture.




Figure 6: User selection in the patient mode of the RPGS

Upon clicking the “Next” button on the bottom right corner of the screen, the patient will be brought to the screen shown in Figure 7. This screen will show a real-time feed from the D435. A guidance skeleton is shown in red initially, while the patient’s joints (called keypoints in Section 7) are highlighted in blue circles. A counter is also present on the top right hand corner showing the number of repetitions completed and the targeted number of repetitions of the current set. Here, the patient will try to attempt the chosen exercise as shown on the user viewer by matching their detected keypoints to that of the guidance skeleton. During the exercise, if the patient has not matched

the posture of the guidance skeleton, or is doing something incorrectly, the guidance skeleton will remain red. Warnings will appear on the screen detailing how the patient should correct their posture. Once the patient completes one repetition of the exercise (consisting of a series of different postures), the completed number of repetitions will increase count by 1.




Figure 7: User viewer in the patient mode of the RPGS

When the user matches the guidance skeleton correctly, the guidance skeleton will flash green and change posture immediately. This is shown in Figure 8.




Figure 8: Guidance skeleton flashing green

4.3.2 Physiotherapist mode

Upon logging into RPGS in physiotherapist mode, the physiotherapist will be brought to the screen shown in Figure 9. The physiotherapist mode consists of two functions, where the physiotherapist is able to view the patient’s medical and application history and also make adjustments to the level of difficulty of the exercise for the patient.




Figure 9: Main page in the physiotherapist mode of the RPGS

After clicking on “Edit exercise”, the physiotherapist will be brought to the screen shown in Figure 10. There, the physiotherapist will be able to customise the exercise difficulty through adjusting the parameters that govern the correctness of each exercise.

One example of how exercises can be customised is to change the angles between relevant keypoints and the respective tolerances of the angles. The number of repetitions per set can also be changed at the bottom left corner.



Figure 10: User editor in the physiotherapist mode of the RPGS

5. Characterisation

Characterisation was carried out after the keypoints were successfully mapped onto a 3D Cartesian coordinate space by deprojection. Characterisation was done to ensure that the experiments were repeatable and the RPGS reliable.

5.1 Deprojection

The 2D coordinates (x, y) of the keypoints extracted from OpenPose will be in pixel coordinates, whereas the depth coordinate (z) will be in meters. This is further elaborated on in Section 7.2. Hence, there is a need to normalize all 3 coordinates to meters for an accurate representation of these points in 3D space. The (x, y) coordinates in pixels have their origin situated at the top left hand corner of the screen, with the positive x -direction towards the right of the screen and the positive y -direction towards the bottom of the screen. This is further complicated by the fact that pixels are wider in real space when they are further from the camera, because the field of view encompasses a larger area. This is shown in Figure 11.

With these problems taken into account, a deprojection algorithm was adopted to normalise the (x, y) pixel coordinates to meters, with the origin of the coordinate system being at the centre of the camera. This will provide us with a representation of the keypoint coordinates (x, y, z) in meters, allowing further characterisation to be done.




Figure 11: Difference between the area covered by the field of view for different distances from the camera¹⁰

5.2 Length of rod

The length of rod characterisation attempts to test the accuracy of measurements done using RPGS in 3D after being able to obtain the (x, y, z) coordinates of each keypoints in meters. The experiment involved the use of a rod that is held by a person at fixed points. The distance between the centre of both wrists were measured using a measuring tape and found to be 0.82m (with an error of $\pm 0.05\text{cm}$). Holding the stick at the same 2 points throughout the entire experiment ensures that the relative distances between the wrist keypoints (4 and 7 according to Figure 20) is invariant (in real space) under the rotation and translation of the rod. Using OpenPose, the keypoints of

both wrists were detected and their depths found using the RealSense SDK.

The keypoints along with their depth values then goes through the deprojection function to obtain the 3D coordinates (x, y, z) in meters. Using equation (1), the distance between two keypoints can be calculated.

$$\text{length } AB = \sqrt{(A_x - B_x)^2 + (A_y - B_y)^2 + (A_z - B_z)^2} \quad \dots (1)$$




Figure 12: Length calculation using a modified version of PRGS. The distance between both wrists are shown on the terminal on the right

As seen in Figure 12, a user can be seen rotating his body along with the stick. Before using the camera system, the actual distance between both his wrists was at 0.82m apart. The user then stood facing the D435 camera and rotated his body to the right, then proceeded to move the rod upwards and downwards. In the terminal on the right, there is a column of values that shows what the camera system and distance algorithm calculated the distance between the wrists keypoints to be.

300 data points were sampled, with the average distance between the 2 keypoints being 0.832m and having a standard deviation of 0.015m, rounded to the nearest 3 decimal places. For purposes of large movements of limbs that physiotherapy is concerned with, this level of accuracy is acceptable and demonstrates that the deprojection algorithm adopted is functional.

5.3 Fluctuations of a fixed point

This experiment involved the characterisation of the stability of the keypoints identified by OpenPose and essentially tests the accuracy of OpenPose in detecting keypoints.

A person was made to sit on a chair with his back against the chair back. This is to ensure that the person remains stationary throughout the course of the experiment. The experimental setup is shown with the person and the camera in Figure 13.



Figure 13: Experimental setup for keypoint fluctuations experiment

The coordinates of the chest keypoint (x, y, z) in meters is recorded every frame and the experiment was terminated after at least 300 frames have elapsed. The standard deviations of the x, y, z points (in meters) are collated in Table 1, rounded off to the nearest 3 decimal places.

| Distance | X std dev | Y std dev | Z std dev |
|----------|-----------|-----------|-----------|
| 1m | 0.006 | 0.004 | 0.093 |
| 2m | 0.006 | 0.004 | 0.093 |
| 3m | 0.006 | 0.004 | 0.092 |

Table 1: standard deviations (in metres) of the x, y and z points at 1m, 2, and 3m respectively

Evidently, the standard deviations for all 3 coordinates do not vary much regardless of distance from the camera. However, the standard deviation of the depth measurement varies slightly. This will not be a problem for the exercise we want to demonstrate in Section 6, since adopting the correct posture will effectively line up all keypoints along a plane and they will not have significant differences in depth. However, further calibration and characterisation would be needed for more complicated postures that involve keypoints being at different depths.

6. Programmed exercise - single arm raise and rotation

For the purpose of this design project, we will be programming a simple exercise that makes use of the 3D imaging function from the D435 and data from the accelerometer. The exercise is shown in Figure 14. Note that for easier visualisation purposes (especially for the wrist area), the wearable module is not worn by the model patient. However, the wearable module is to be worn by the patient throughout the entire course of the exercise to detect for wrist rotation.




Figure 14: Pre-programmed exercise: Single Arm Raise and Rotation

The exercise involves 4 steps:

1. The patient first adopts the starting posture by having their hands by their sides, with their palms facing their thighs.
2. The patient raises their right arm such that it is 90° to their body, and 90° to the direction of view of the camera. The patient's palm should be facing the floor
3. The patient rotates their right arm along its own axis by 180° while maintaining the same posture as step 2. The patient's palm should be facing the ceiling.
4. The patient returns to the posture in step 1. This counts as 1 repetition of the exercise.

Throughout this exercise, various conditions have been implemented to detect for error in posture. Some examples include (but are not limited to): the palms being oriented the wrong way in steps 2 and 3, and the arm not being raised to 90° in step 2. A full list of these conditions, possible errors, and how they are implemented within the software will be detailed in Section 7.3.

7. Core Architecture

In this section, the underlying core architecture of the software component of the RPGS will be explained. The core architecture is made up of two toolkits: the RealSense SDK¹¹ (Software Development Kit) and OpenPose¹². The RealSense SDK is a Github real-time multi-person keypoint detection library for body, face, hands, and foot estimation. It is only able to detect points on a 2D flat plane with no depth elements involved. The mapping of 2D points on a plane to 3D space using deprojection was already explained in Section 5.1

A pictorial representation of the software architecture is shown in Figure 15 below.




Figure 15: An overview of the technical workflow of RPGS

As seen in Figure 15, after the patient selects the appropriate exercise to do, the selected python script will first initialise both RealSense and OpenPose, then execute the RealSense SDK's pipeline, before running OpenPose within the SDK's pipeline. Conditions to decide whether the patient's posture is correct or not are also verified within the pipeline after extracting the necessary

information from the D435 camera, RealSense SDK and OpenPose. The pipeline including OpenPose then runs continuously until the user chooses to exit the exercise. The following sections will discuss the technical aspects of the RealSense SDK and OpenPose and how the correctness of poses are calculated.

7.1 RealSense SDK

The RealSense SDK works synergistically with the D435 camera to allow for image and depth manipulation. The native coding language of the SDK is in C++, but it also includes wrappers for other languages like Python, which was chosen to be the language to write RPGS. A pictorial representation of the SDK's pipeline can be seen in Figure 16 below.




Figure 16: Breakdown of the operations happening within the RealSense SDK pipeline

As seen in Figure 16, the RealSense SDK will first undergo initialisation before opening up the pipeline. The initialisation will first enable both colour and depth sensors to collect a colour and depth stream respectively at a resolution of 640x360. The aforementioned resolution was chosen as it is the most common aspect ratio¹³, and the experimental setup uses a laptop with the same aspect ratio. A resolution of 360p was chosen to improve the frames per second that is able to be shown on the screen. After setting up the sensors and resolution, the last step involves creating a function that will be used later during step 2 to align the depth frame to the colour frame. This is because the depth frame has a larger field of view than the colour frame.



Figure 17: Field of view comparison between the colour (left) and depth (right) frame

As seen in Figure 17, the depth frame has a larger field of view when compared to the colour image. This is evident in the head of the person being imaged in the depth frame, but not in the colour frame. Hence, there is a need to make sure both frames contain the same field of view as accurate depth calculation would require the frames to be aligned

When the pipeline starts running, the SDK will first extract a frame from the D435 camera stream and split it into a colour and depth frame for further processing. After aligning both frames, the SDK then extracts the colour frame's intrinsic information, such as the dimensions of the camera image, principal point, focal points, distortion model and distortion coefficients¹⁴. This is because the intrinsic information of the colour frame would be necessary to create the final 3D coordinates of any pixel point. In step 4, the 2D colour image will be saved for OpenPose to process. After extracting the 2D coordinates of the patient's joint coordinates, also known as keypoints, calculations would be made (based on the customisation) from those keypoints to determine whether the patient's posture is correct or not. In Step 7, information from the wearable module would be sampled before drawing the user skeleton, guidance skeleton, any applicable warnings and a help function. Finally, as long as the user does not close or logout from RPGS, the pipeline will continuously loop from step one to ten.

7.2 OpenPose

OpenPose is an open source, real-time keypoint detection library, which is able to detect the human body, face, hand and feet. It is able to detect keypoints either in 2D, using a single camera or in 3D, using multiple cameras. For the user's convenience of not having to set up more than one camera and budget constraints, RPGS will first use OpenPose to obtain the 2D keypoints of the user's body, before subsequently lifting the keypoints to 3D by using the RealSense SDK. This is detailed in Section 7.3. Hence,

OpenPose was written in Python as well for language compatibility. An illustration of OpenPose's workflow can be seen in Figure 18 below.




Figure 18: Breakdown of the operations happening when running OpenPose

In Figure 18, OpenPose will also first undergo an initialisation stage. It will first load the chosen pose model called “body_25”, as it is faster and more accurate than the other two options called “COCO” and “MPI”¹⁵. An illustration of the body_25 model is shown in Figure 19 below.




Figure 19: Labelled keypoint of body_25 model¹⁶. Image is mirrored

In Figure 19, the body_25 model involves identifying 25 different keypoints around the human body, and theoretically, we can program different exercises involving any of these 25 different keypoints. However, since the chosen exercise in Section 6 only involves the arms, RPGS uses a modified version of body_25 as a representation for the user and guidance skeleton as seen in Figure 20 below

(the model used to identify the keypoints is still the full version of body_25). However, it is still possible to return the keypoints for the ankles and face if they are needed to be programmed for more exercises focusing on those keypoints.




Figure 20: Modified version of body_25 used in RPGS

As seen in Figure 20, the modified version removes ten redundant keypoints from the head and from the feet. This makes it look simpler and will be easier to understand.

Another important parameter to initialise is the number of people OpenPose will take into account. Since RPGS is meant for one user at a time, OpenPose will only process the most prominent person it can detect. This parameter could incur better performance since computation is only done on one person, and any false positive detection in the background is eliminated. Afterwhich, OpenPose is set to asynchronous mode, as it will be processing an image and not the default camera stream which is occupied by the RealSense SDK. In the pipeline, OpenPose will save the computed keypoints into a variable in the code, which is more efficient than the previous workflow as seen in Figure 21 below.




Figure 21: Previous version of the workflow of OpenPose

As seen in Figure 21, the previous version of RPGS ran OpenPose in synchronous mode and had to save the computed keypoints into a json file. This caused RPGS to run two additional functions to extract the relevant 2D keypoints from OpenPose. Hence, the current implementation of OpenPose is much more efficient than before.

7.3 Calculation of Conditions

Before deciding whether the user's pose is considered correct or not, calculations need to be based on the pre-programmed conditions and whether the relative orientations of the detected keypoints fit those conditions. For the purposes of our chosen exercise, the pre-programmed conditions will only involve the angles between various keypoints, and the spatial orientation of the wearable module.

7.3.1 Calculation of Angles

For the chosen exercise, the keypoints identified to be relevant are the:

1. Chest
2. Right shoulder
3. Right elbow
4. Right wrist
8. Centre of waist
9. Right pelvis
12. Left pelvis

Figure 22 shows a pictorial representation of the relevant keypoints and the angles. For example, angle 1-2-3 simply refers to the smaller angle (less than 180°) between the keypoints 1, 2, 3. Vector N is defined as the normal of the plane defined by the keypoints 1, 9, 12. It is chosen to be

oriented out of the page, or towards the front of the patient, by taking the cross product of angles 1-9 and 1-12. Vector 1-3 refers to the vector in the originating from keypoint 1 and ending at keypoint 3. All relevant vectors have been normalized to unit vectors before calculating angles, and a standard tolerance of $\pm 10^\circ$ is allocated to account for the differences in the skeletal structure of patients and possible systematic errors such as the points 1, 2, 3, 4 not lying exactly on a plane because keypoint 1 (chest) will be closer to the camera than keypoints 2 (shoulder), 3 (elbow), 4 (wrist) due to its protruding nature. However, the physiotherapist in physiotherapist mode can alter the ideal angle and respective tolerances as well (Figure 10).

Throughout the exercise, we pre-programmed five conditions, being:

1. Angle 1-2-3 = $(180 \pm 10)^\circ$
2. Angle 2-3-4 = $(180 \pm 10)^\circ$
3. Angle 8-1-4 = $(90 \pm 10)^\circ$
4. Angle between vector N and vector 1-3 = $(90 \pm 10)^\circ$
5. Angle between vector N and vector 1-4 = $(90 \pm 10)^\circ$




Figure 22: Pictorial representation of keypoints and angles

The calculation of these angles is achieved by first finding the depth of keypoints 1, 3, 4, 8, 9 and 12 in meters using the built-in function in the RealSense SDK. Next, all the keypoints pass through another built-in function called deprojection. According to Intel's documentation, deprojection takes a 2D coordinate on a stream's image, as well as its depth, and maps it to a 3D point location within the stream's associated 3D coordinate space¹⁷. This basically converts the 2D coordinates from OpenPose into 3D coordinates in the camera's coordinate space. Once the 3D coordinates are found, the vectors 1-3, 1-4, 1-9 and 1-12 can then be computed along with their lengths using

$$\text{Vector } AB = (B_x - A_x, B_y - A_y, B_z - A_z) \quad \dots (2)$$

$$\text{length } AB = \sqrt{(A_x - B_x)^2 + (A_y - B_y)^2 + (A_z - B_z)^2} \quad \text{--- (3)}$$

Where $A_{x/y/z}$ is keypoint A in the x, y and depth axis, and $B_{x/y/z}$ is keypoint B in the x, y and depth axis. Finding both values allows for the computation of a normalised vector by dividing the vector with its length. Next, the normal vector is calculated using the cross product function provided by the numpy module between normalised vectors 1-9 and 1-12. The normal vector would also be normalised. Finally, the angles N-1-3 and N-1-4 can be found by using the dot product function provided in the numpy module, between the normal vector and the normalised vectors 1-3 and 1-4 respectively.

The angles 1-2-3, 2-3-4 and 8-1-4 are 2D angles also calculated using geometry methods. In this case, the vectors and lengths can be calculated using

$$\text{Vector } AB = (B_x - A_x), (B_y - A_y) \quad \text{--- (4)}$$

$$\text{length } AB = \sqrt{(A_x - B_x)^2 + (A_y - B_y)^2} \quad \text{--- (5)}$$

For example, the angle 1-2-3 can be found by dividing equation (3) and (4) to find the normalised vector 2-1 and 2-3. Then angle 1-2-3 can be found by using the dot product function between normalised vectors 2-1 and 2-3.

If any of the keypoints mentioned above is not detected, the angles that contain that keypoints would equate to zero.

7.3.2 Accelerometer conditions

Data readings from the accelerometer in the wearable module is also used as a condition in RPGS. A copy of the Arduino code can be found in Appendix III.

Float variables (accelX, accelY, accelZ) for the acceleration in x, y, z directions according to the coordinate system defined in Figure 4 are first initialised. A loop is then activated for the accelerometer registers to start recording data. The data is transferred to the Arduino UNO with the Wire.h library and normalized to the acceleration due to gravity (9.81ms^{-2}). A serial.print command is activated for each of the float variables to be displayed on the serial monitor. For purposes of our chosen exercise, data from the wearable module is only used to determine the spatial orientation of itself. That is to say, most of the acceleration the wearable module experiences is due to gravity.

After uploading the Arduino code into the Arduino, the acceleration data can be imported directly into the RealSense pipeline using the pyserial module. The data is then filtered to identify the highest absolute value out of all 3 orthogonal directions, and the axis that corresponds to that value. This is because the wearable module will be most aligned to the axis that experiences the greatest magnitude of gravitational acceleration. With that, RPGS will be able to detect the spatial orientation of the wearable module and thus, be able to use it as another condition to check if the patient is rotating their arm to the correct orientation.

The accelerometer conditions for all 4 steps are:

1. Predominately negative y-direction
2. Predominately negative x-direction
3. Predominately positive x-direction
4. Predominately negative x-direction (same as 1)

7.3.3 Deciding the outcome

After checking against the conditions mentioned previously, the correctness of the posture will be shown through the colour of the guidance skeleton. When the guidance skeleton flashes green, it means the pose is correct, and the guidance skeleton will change posture to the next step of the exercise. However, if the guidance skeleton turns red, it means the pose is incorrect. An illustration of the sequence of checks can be seen in Figure 23 below.




Figure 23: Workflow of how the outcome of the user's pose is decided using the guidance skeleton

As seen in Figure 23, after a physiotherapist sets the angles and the respective tolerances, the values will be saved into a json file. When the RealSense SDK's pipeline reaches step 9, it will compare the calculated values with their

corresponding saved values from those json files. If the calculated angles fall within the range of values saved and the rotational direction is correct for that pose, the guidance skeleton will be drawn in green to indicate that the pose is correct. On the other hand, if any of the conditions fails, the guidance skeleton will remain as red to indicate that the pose is incorrect. A warning system is also in place in the event that the outcome is incorrect to guide the user to the correct pose.

A list of warnings for the Single Arm Raise and Rotation exercise is shown below:

Distance from camera:

1. “Move forward/backward” when any of the keypoints are <2m or >3m

Arm movement (steps 2 and 3):

1. “Move your arm up” when angle 8-1-4 is less than the lower bound of the tolerance
2. “Move your arm down” when angle 8-1-4 is more than the upper bound of the tolerance
3. “Move your arm backwards” when either the angle between vector N and vector 1-3 or vector 1-4 is less than the lower bound of the tolerance
4. “Move your arm forwards” when either the angle between vector N and vector 1-3 or vector 1-4 is less than the upper bound of the tolerance
5. “Straighten your elbow” when angle 2-3-4 is out of the range of the tolerance

Rotation (all steps):

1. “Rotate your hand until your palm faces your leg” (step 1 & 4) when the wearable module is not oriented sufficiently in the negative y-direction
2. “Rotate your hand until your palm faces the floor” when the wearable module is not oriented sufficiently in the negative x-direction
3. “Rotate your hand until your palm faces the ceiling” when the wearable module is not oriented sufficiently in the positive x-direction

Lastly, by leveraging on the condition checking procedures, RPGS is also able to count the number of repetitions the patient has completed, against the total number of repetitions they are to do, as set by the physiotherapist.

8. Results and Discussion

The RPGS was able to demonstrate static posture recognition with the Single Arm Raise and Rotation in front of a live ZOOM audience during the ESP3902 Final Presentation. 3 different programs were run: purely 2D posture recognition, purely 3D posture recognition, and 3D posture recognition with the wearable module. This allowed the illustration of how the purely 2D and 3D posture recognition programs were unable to recognise some errors in posture due to their limited degrees of freedom. This included the inability to define a normal vector N in 2D, and the inability to detect rotation of joints in 3D. These problems were solved and demonstrated by using the 3D posture recognition with the wearable module.

However, feedback given included the response time for detecting the correct posture being relatively slow. This has to do with the Arduino within the wearable module transferring through PySerial at a faster rate than the frames per second, causing a misalignment between the frame and the specific values of the acceleration corresponding to that frame. Improvements can be made in the form of not sampling every single frame and introducing a common clocking system to sync both data from the Arduino and the frame rate.

All in all, the RPGS was able to execute a proof-of-concept exercise by using the functions of the 3D camera and accelerometer. This shows that more complicated exercises can be programmed in future, using the same fundamental concepts that have been explained in earlier sections.

9. Future Works

Although the RPGS is currently able to function for a simple exercise of a single arm raise followed by rotation of the arm, there are still further improvements that can be made in future with more time.

Firstly, the frame rate can be increased from 8 to 15 which is the frame rate of Assistive Rehab¹⁸ which is another physiotherapy guidance platform that involves the patient following the actions of a life-sized robot. Assistive Rehab uses both OpenPose and RealSense libraries together with the physical robot to guide the patient on how to attempt the exercise correctly. This will make RPGS even more accurate in posture detection and have a faster response time.

Secondly, dynamic monitoring of velocity of the limbs can also be implemented on top of having static poses. This will give both the patient an accurate analysis of the speed of their limb movements, effectively adding an additional degree of freedom to quantify and analyse. Possible applications of this can be in the form of detecting how fast the exercise is being done, to determine if the speed of the exercise is too fast or too slow.

Lastly, machine learning can possibly be implemented for complicated exercises poses that can be difficult for conditions to be set. This will involve multiple repetitions of the exercise by different patients to accumulate a diverse and large data set consisting of time series data of the joints. Time series analysis can then be done in real time to identify if the patient is performing the exercise accurately.

10. Conclusion

The RPGS is a home-based physiotherapy guidance system that is easy to set up and operate. It has proven to be accurate in determining the correctness of exercises, and able to provide specific real-time feedback of how to correct one's pose when one is performing it incorrectly. Through the proof-of-concept demonstration of the Single Arm Raise and Rotation exercise, the RPGS has been shown to be able to leverage on the functionalities of both D435 3D camera and accelerometer to detect the correctness of postures based on the spatial relationships between keypoints in 3D and spatial orientation of the wearable module.

11. References

- [1] What is physiotherapy? (n.d.). Retrieved November 15, 2020, from <https://www.csp.org.uk/careers-jobs/what-physiotherapy>
- [2] The Importance of Doing Physio Exercises. (2020, April 06). Retrieved November 15, 2020, from <https://www.cairnhillphysiotherapy.co.nz/blog/the-importance-of-doing-physio-exercises/>
- [3] Tang, R., Yang, X., Bateman, S., Jorge, J., & Tang, A. (2015). Physio@Home. *Proceedings of the 33rd Annual ACM Conference on Human Factors in Computing Systems - CHI '15*. doi:10.1145/2702123.2702401
- [4] Graves, J. M., Iyer, K. R., Willis, M. M., Ebel, B. E., Rivara, F. P., & Vavilala, M. S. (2013). Emergency department-reported injuries associated with mechanical home exercise equipment in the USA. *Injury Prevention*, 20(4), 281-285. doi:10.1136/injuryprev-2013-040833
- [5] Geskey, J. M. (2011). Home Exercise Equipment-Related Injuries. *AAP Grand Rounds*, 26(3), 28-28. doi:10.1542/gr.26-3-28
- [6] Depth Camera D435. (2020, September 25). Retrieved November 15, 2020, from <https://www.intelrealsense.com/depth-camera-d435/>
- [7] Arduino Uno Rev3. (n.d.). Retrieved November 15, 2020, from <https://store.arduino.cc/usa/arduino-uno-rev3>
- [8] MPU-6050: TDK. (n.d.). Retrieved November 15, 2020, from <https://invensense.tdk.com/products/motion-tracking/6-axis/mpu-6050/>
- [9] What is PyQt? (n.d.). Retrieved November 15, 2020, from <https://riverbankcomputing.com/software/pyqt/intro>
- [10] Pourreza Ghoushchi, Vahid. (2015). Opto-Mechanical Design And Development of an Optodigital Confocal Microscope. 10.13140/RG.2.2.33261.38884.
- [11] Developers area – Intel RealSense SDK 2.0 - Depth and Tracking cameras. (2020, October 22). Retrieved November 15, 2020, from <https://www.intelrealsense.com/developers/>
- [12] CMU-Perceptual-Computing-Lab. (n.d.). CMU-Perceptual-Computing-Lab/openpose. Retrieved November 15, 2020, from <https://github.com/CMU-Perceptual-Computing-Lab/openpose>
- [13] Rusen, C. A. (2020, October 25). What do the 720p, 1080p, 1440p, 2K, 4K resolutions mean? What are the aspect ratio & orientation? Retrieved November 15, 2020, from <https://www.digitalcitizen.life/what-screen-resolution-or-aspect-ratio-what-do-720p-1080i-1080p-mean>

- [14] Pterneas, V. (2018, November 02). Intel RealSense D415/435: Coordinate Mapping in C#. Retrieved November 15, 2020, from
<https://lightbuzz.com/intel-realsense-coordinate-mapping/>
- [15] CMU-Perceptual-Computing-Lab. (2020, November 09). CMU-Perceptual-Computing-Lab/openpose. Retrieved November 15, 2020, from
<https://github.com/CMU-Perceptual-Computing-Lab/openpose/blob/master/doc/faq.md>
- [16] CMU-Perceptual-Computing-Lab. (2019, November 10). CMU-Perceptual-Computing-Lab/openpose. Retrieved November 15, 2020, from
<https://github.com/CMU-Perceptual-Computing-Lab/openpose/blob/master/doc/output.md>
- [17] Projection in Intel RealSense SDK 2.0. (2019). Retrieved November 15, 2020, from
<https://dev.intelrealsense.com/docs/projection-in-intel-realsense-sdk-20>
- [18] Vasco, V. (n.d.). Assistive-rehab. Retrieved November 16, 2020, from
<https://robotology.github.io/assistive-rehab/doc/mkdocs/site/>

Appendix I

(GUI code)

(main UI code)

```
1 #####  
2 ##  
3 ## BY: WANDERSON M.PIMENTA  
4 ## PROJECT MADE WITH: Qt Designer and PySide2  
5 ## V: 1.0.0  
6 ##  
7 #####  
8  
9 import sys  
10 import platform  
11 from PyQt5 import QtCore, QtGui, QtWidgets  
12 from PyQt5.QtCore import (QCoreApplication, QPropertyAnimation, QDate, QDateTime, QMetaObject, QObject, QPoint, QRect, QSize, QTime, QUrl, Qt, QEvent)  
13 from PyQt5.QtGui import (QBrush, QColor, QConicalGradient, QCursor, QFont, QFontDatabase, QIcon, QKeySequence, QLinearGradient, QPalette, QPainter, QPixmap, QRadialGradient)  
14 from PyQt5.QtWidgets import *  
15  
16 ## ==> SPLASH SCREEN  
17 from ui_splash_screen import Ui_SplashScreen  
18  
19 ## ==> GLOBALS  
20 counter = 0  
21  
22  
23  
24  
25 # SPLASH SCREEN  
26 class SplashScreen(QMainWindow):  
27     def __init__(self):  
28         QMainWindow.__init__(self)  
29         self.ui = Ui_SplashScreen()  
30         self.ui.setupUi(self)  
31  
32         ## UI ==> INTERFACE CODES  
33         #####  
34  
35         ## REMOVE TITLE BAR  
36         self.setWindowFlag(QtCore.Qt.FramelessWindowHint)  
37         self.setAttribute(QtCore.Qt.WA_TranslucentBackground)  
38  
39  
40         ## DROP SHADOW EFFECT  
41         self.shadow = QGraphicsDropShadowEffect(self)  
42         self.shadow.setBlurRadius(20)  
43         self.shadow.setXOffset(0)
```

```

44     self.shadow.setYOffset(0)
45     self.shadow.setColor(QColor(0, 0, 0, 60))
46     self.ui.dropShadowFrame.setGraphicsEffect(self.shadow)
47
48     ## QTIMER ==> START
49     self.timer = QtCore.QTimer()
50     self.timer.timeout.connect(self.progress)
51     # TIMER IN MILLISECONDS
52     self.timer.start(35)
53
54     # CHANGE DESCRIPTION
55
56     # Initial Text
57     self.ui.label_description.setText("<strong>WELCOME</strong> TO RPGS!")
58
59     # Change Texts
60     QtCore.QTimer.singleShot(3500, lambda: self.ui.label_description.setText("<strong>LOADING</strong> DATABASE"))
61     QtCore.QTimer.singleShot(5000, lambda: self.ui.label_description.setText("<strong>LOADING</strong> USER INTERFACE"))
62     QtCore.QTimer.singleShot(6500, lambda: self.ui.label_description.setText("<strong>LOADING</strong> Scripts"))
63
64
65     ## SHOW ==> MAIN WINDOW
66     #####
67     self.show()
68     ## ==> END ##
69
70 ## ==> APP FUNCTIONS
71 #####
72 def progress(self):
73
74     global counter
75
76     # SET VALUE TO PROGRESS BAR
77     self.ui.progressBar.setValue(counter)
78
79     # CLOSE SPLASH SCREE AND OPEN APP
80     if counter > 100:
81         # STOP TIMER
82         self.timer.stop()
83
84         # SHOW MAIN WINDOW
85         import RPGS
86
87         # CLOSE SPLASH SCREEN
88         self.close()
89
90         # INCREASE COUNTER
91         counter += 0.65
92
93
94
95     if __name__ == "__main__":
96         app = QApplication(sys.argv)
97         window = SplashScreen()
98         sys.exit(app.exec_())
99
100

```

(loading bar code)

```
1  # -*- coding: utf-8 -*-
2
3 ##### Form generated from reading UI file 'splash_screenXBSmkq.ui'
4 ## Form generated from reading UI file 'splash_screenXBSmkq.ui'
5 ##
6 ## Created by: Qt User Interface Compiler version 5.14.1
7 ##
8 ## WARNING! All changes made in this file will be lost when recompiling UI file!
9 #####
10
11 from PyQt5.QtCore import (QCoreApplication, QMetaObject, QObject, QPoint,
12     QRect, QSize, QUrl, Qt)
13 from PyQt5.QtGui import (QBrush, QColor, QConicalGradient, QCursor, QFont,
14     QFontDatabase, QIcon, QLinearGradient, QPalette, QPainter, QPixmap,
15     QRadialGradient)
16 from PyQt5.QtWidgets import *
17
18
19 class Ui_SplashScreen(object):
20     def setupUi(self, SplashScreen):
21         if SplashScreen.objectName():
22             SplashScreen.setObjectName(u"SplashScreen")
23             SplashScreen.resize(680, 400)
24             self.centralwidget = QWidget(SplashScreen)
25             self.centralwidget.setObjectName(u"centralwidget")
26             self.verticalLayout = QVBoxLayout(self.centralwidget)
27             self.verticalLayout.setSpacing(0)
28             self.verticalLayout.setObjectName(u"verticalLayout")
29             self.verticalLayout.setContentsMargins(10, 10, 10, 10)
30             self.dropShadowFrame = QFrame(self.centralwidget)
31             self.dropShadowFrame.setObjectName(u"dropShadowFrame")
32             self.dropShadowFrame.setStyleSheet(u"QFrame { \n"
33             "    background-color: rgb(56, 58, 89); \n"
34             "    color: rgb(220, 220, 220); \n"
35             "    border-radius: 10px; \n"
36             "}")
37             self.dropShadowFrame.setFrameShape(QFrame.StyledPanel)
38             self.dropShadowFrame.setFrameShadow(QFrame.Raised)
39             self.label_title = QLabel(self.dropShadowFrame)
40             self.label_title.setObjectName(u"label_title")
41             self.label_title.setGeometry(QRect(0, 90, 661, 61))
42             font = QFont()
43             font.setFamily(u"Segoe UI")
```

```

44     font.setPointSize(20)
45     self.label_title.setFont(font)
46     self.label_title.setStyleSheet(u"color: rgb(254, 121, 199);")
47     self.label_title.setAlignment(Qt.AlignCenter)
48     self.label_description = QLabel(self.dropShadowFrame)
49     self.label_description.setObjectName(u"label_description")
50     self.label_description.setGeometry(QRect(0, 150, 661, 31))
51     font1 = QFont()
52     font1.setFamily(u"Segoe UI")
53     font1.setPointSize(14)
54     self.label_description.setFont(font1)
55     self.label_description.setStyleSheet(u"color: rgb(98, 114, 164);")
56     self.label_description.setAlignment(Qt.AlignCenter)
57     self.progressBar = QProgressBar(self.dropShadowFrame)
58     self.progressBar.setObjectName(u"progressBar")
59     self.progressBar.setGeometry(QRect(50, 280, 561, 23))
60     self.progressBar.setStyleSheet(u"QProgressBar {\n"
61     " \n"
62     " background-color: rgb(98, 114, 164);\n"
63     " color: rgb(200, 200, 200);\n"
64     " border-style: none;\n"
65     " border-radius: 10px;\n"
66     " text-align: center;\n"
67     "}\n"
68     "QProgressBar::chunk{\n"
69     " border-radius: 10px;\n"
70     " background-color: qlineargradient(spread:pad, x1:0, y1:0.511364, x2:1, y2:0.523, stop:0 rgba(254, 121, 199, 255), stop:1 rgba(170, 85, 255, 255));\n"
71     "}")
72     self.progressBar.setValue(24)
73     self.label_loading = QLabel(self.dropShadowFrame)
74     self.label_loading.setObjectName(u"label_loading")
75     self.label_loading.setGeometry(QRect(0, 320, 661, 21))
76     font2 = QFont()
77     font2.setFamily(u"Segoe UI")
78     font2.setPointSize(12)
79     self.label_loading.setFont(font2)
80     self.label_loading.setStyleSheet(u"color: rgb(98, 114, 164);")
81     self.label_loading.setAlignment(Qt.AlignCenter)
82     self.label_credits = QLabel(self.dropShadowFrame)
83     self.label_credits.setObjectName(u"label_credits")
84     self.label_credits.setGeometry(QRect(20, 350, 621, 21))
85     font3 = QFont()
86     font3.setFamily(u"Segoe UI")
87
88     font3.setPointSize(10)
89     self.label_credits.setFont(font3)
90     self.label_credits.setStyleSheet(u"color: rgb(98, 114, 164);")
91     self.label_credits.setAlignment(Qt.AlignRight|Qt.AlignTrailing|Qt.AlignVCenter)
92
93     self.verticalLayout.addWidget(self.dropShadowFrame)
94
95     SplashScreen.setCentralWidget(self.centralwidget)
96
97     self.retranslateUi(SplashScreen)
98
99     QMetaObject.connectSlotsByName(SplashScreen)
100
101     # setupUi
102
103     def retranslateUi(self, SplashScreen):
104         SplashScreen.setWindowTitle(QCoreApplication.translate("SplashScreen", u"MainWindow", None))
105         self.label_title.setText(QCoreApplication.translate("SplashScreen", u" RealSense Physiotherapy Guidance System", None))
106         self.label_description.setText(QCoreApplication.translate("SplashScreen", u"<strong>APP</strong> DESCRIPTION", None))
107         self.label_loading.setText(QCoreApplication.translate("SplashScreen", u"loading...", None))
108
109         # retranslateUi

```

(RPGS UI code)

```
 1 import sys
 2 import os
 3 import json
 4
 5 from PyQt5 import uic, QtWidgets, QtCore
 6 from PyQt5.QtWidgets import *
 7 from PyQt5.QtCore import *
 8 from PyQt5.QtGui import *
 9 import webbrowser # to open pdf files
10 from PyQt5.QtGui import QMovie, QPixmap # for the gif
11
12
13 qtCreatorFile = "Main_GUI.ui"
14
15
16 class Ui(QtWidgets.QMainWindow):
17
18     def __init__(self):
19         super(Ui, self).__init__()
20         uic.loadUi(qtCreatorFile, self)
21
22         # REMOVE TITLE BAR
23         self.setWindowFlag(QtCore.Qt.FramelessWindowHint)
24         self.setAttribute(QtCore.Qt.WA_TranslucentBackground)
25
26     def moveWindow(self, event):
27         # IF LEFT CLICK MOVE WINDOW
28         if event.buttons() == Qt.LeftButton:
29             self.move(self.pos() + event.globalPos() - self.dragPos)
30             self.dragPos = event.globalPos()
31             event.accept()
32
33         # SET TITLE BAR
34         self.title_bar.mouseMoveEvent = moveWindow
35
36         # CLOSE
37         self.btn_close.clicked.connect(lambda: self.close())
38
39         # MINIMIZE
40         self.btn_minimize.clicked.connect(lambda: self.showMinimized())
41
42         # change the program's main window name
43         self.setWindowTitle('RPGS')
```

```

44     # Just for [3] label on top left
45     self.label_text = ""
46
47     # Prevent resizing
48     self.setFixedSize(self.size())
49
50
51
52     # Login page
53     self.stackedWidget.setCurrentIndex(0)
54     self.Username_field.returnPressed.connect(self.Login_btn.click) # enter key act as clicking the login button
55     self.Password_field.returnPressed.connect(self.Login_btn.click)
56     self.Login_btn.clicked.connect(self.loginfunction) # main login function
57
58     # Choose exercise page (PATIENT) [1]
59     self.Exercise_ChooseExercise.addItem("Half Arm Raise [2D]") # add the items to be selectable
60     self.Exercise_ChooseExercise.addItem("Half Arm Raise [3D]")
61     self.Exercise_ChooseExercise.addItem("Half Arm Raise [3D + Arduino]")
62     self.Exercise_ChooseExercise.itemClicked.connect(self.Exercise_Display) # run video when exercise selected
63     self.Exercise_NextBtn.clicked.connect(self.run_exercise) # run the py script depending on what user chose
64     self.Logout_btn.clicked.connect(self.logout)
65
66     # Choose exercise page (DEVELOPER) [2]
67     self.Exercise_ChooseExercise_2.addItem("Half Arm Raise [2D]") # add the items to be selectable
68     self.Exercise_ChooseExercise_2.addItem("Half Arm Raise [3D]")
69     self.Exercise_ChooseExercise_2.addItem("Half Arm Raise [3D + Arduino]")
70     self.logout_btn_2.clicked.connect(self.logout)
71     self.edit_exercise_btn.clicked.connect(self.click_to_label_variable) # go to next page while changing label text
72
73     # Dev chosen exercise page [3]
74     self.stackedWidget_edit.setCurrentIndex(0)
75
76     with open(r"D:\STUDIES\YEAR 3 - SEM 1\ESP3902\openpose-master\build\examples\tutorial_api_python\Dev\Json\start.json") as start_json:
77         start = json.load(start_json)
78         for start_123 in start['start_123']:
79             self.start_123_angle.setText(str(start_123["angle"]))
80             self.start_123_upper.setText(str(start_123["upper"]))
81             self.start_123_lower.setText(str(start_123["lower"]))
82         for start_234 in start['start_234']:
83             self.start_234_angle.setText(str(start_234["angle"]))
84             self.start_234_upper.setText(str(start_234["upper"]))
85             self.start_234_lower.setText(str(start_234["lower"]))
86

```

```

87     with open(r"D:\STUDIES\YEAR 3 - SEM 1\ESP3902\openpose-master\build\examples\tutorial_api_python\Dev_Json\end.json") as end_json:
88         end = json.load(end_json)
89         for end_234 in end['end_234']:
90             self.end_234_angle.setText(str(end_234["angle"]))
91             self.end_234_upper.setText(str(end_234["upper"]))
92             self.end_234_lower.setText(str(end_234["lower"]))
93         for end_814 in end['end_814']:
94             self.end_814_angle.setText(str(end_814["angle"]))
95             self.end_814_upper.setText(str(end_814["upper"]))
96             self.end_814_lower.setText(str(end_814["lower"]))
97         for end_13 in end['end_13']:
98             self.end_13_angle.setText(str(end_13["angle"]))
99             self.end_13_upper.setText(str(end_13["upper"]))
100            self.end_13_lower.setText(str(end_13["lower"]))
101        for end_14 in end['end_14']:
102            self.end_14_angle.setText(str(end_14["angle"]))
103            self.end_14_upper.setText(str(end_14["upper"]))
104            self.end_14_lower.setText(str(end_14["lower"]))
105
106    with open(r"D:\STUDIES\YEAR 3 - SEM 1\ESP3902\openpose-master\build\examples\tutorial_api_python\Dev_Json\repetition.json") as rep_json:
107        rep = json.load(rep_json)
108        self.Number_of_reps.setValue(int(rep['repetition'][0]['repetition_number']))
109
110
111    self.start_pose_edit.clicked.connect(lambda: self.stackedWidget_edit.setCurrentIndex(0))
112    self.end_pose_edit.clicked.connect(lambda: self.stackedWidget_edit.setCurrentIndex(1))
113    self.start_save_btn.clicked.connect(self.save_start)
114    self.end_save_btn.clicked.connect(self.save_end)
115    self.dev_edit_back_btn.clicked.connect(lambda: self.stackedWidget.setCurrentIndex(2))
116
117
118    # For the menu bar
119    self.btn_help.clicked.connect(self.Help)
120
121
122    # DONT TOUCH
123    self.show()
124
125 #-----
126
127    # FUNCTIONS FROM HERE ONWARDS
128    def loginfunction(self):
129        username = self.Username_field.text()

```

```

130     password = self.Password_field.text()
131     self.Username_field.clear()
132     self.Password_field.clear()
133     if(username == "patient" and password == "patient"):
134         self.stackedWidget.setCurrentIndex(1)
135     elif(username== "dev" and password == "dev"):
136         self.stackedWidget.setCurrentIndex(2)
137     elif(username == "" and password == ""):
138         return
139     else:
140         error_screen = QMessageBox()
141         error_screen.setWindowTitle("Wrong Login")
142         error_screen.setIcon(QMessageBox.Warning)
143         error_screen.setText("Wrong Username or Password")
144         x = error_screen.exec_()
145
146
147     def Exercise_Display(self):
148         if(self.Exercise_ChooseExercise.currentRow() == 0 or 1 or 2):
149             movie = QMovie('instructional-video-raw.gif')
150             self.Exercise_DisplayExercise.setMovie(movie)
151             movie.start()
152
153
154     def run_exercise(self):
155         listItems = self.Exercise_ChooseExercise.selectedItems()
156         if not listItems:
157             return
158         else:
159             exercise_chosen = self.Exercise_ChooseExercise.currentItem().text()
160             if(exercise_chosen == "Half Arm Raise [2D]"):
161                 exec(open("realsense_OP_Half_arm_raise_2D.py").read())
162             elif(exercise_chosen == "Half Arm Raise [3D]"):
163                 exec(open("realsense_OP_Half_arm_raise_3D_V2.py").read())
164             elif(exercise_chosen == "Half Arm Raise [3D + Arduino]"):
165                 exec(open("realsense_OP_Half_arm_raise_3D_Arduino_V4(for GUI).py").read())
166
167
168     def logout(self):
169         logout_screen = QMessageBox()
170         user_click = logout_screen.question(self,'Log Out', "Are You Sure You Want To Log Out?", logout_screen.Yes | logout_screen.No)
171         if(user_click == logout_screen.Yes):
172             self.stackedWidget.setCurrentIndex(0)

```

```

173
174
175     def Help(self):
176         # CHANGE THE PATH ACCORDING TO THE CORRECT PDF
177         webbrowser.open_new('instruction_manual.pdf')
178
179
180     def click_to_label_variable(self):
181         listItems = self.Exercise_ChooseExercise_2.selectedItems()
182         if not listItems:
183             return
184         else:
185             self.label_text = self.Exercise_ChooseExercise_2.currentItem().text()
186             self.stackedWidget.setCurrentIndex(3)
187             self.Dev_exercise_label.setText(self.label_text)
188
189
190     def save_start(self):
191         save = QMessageBox()
192         user_click = save.question(self,'Starting Pose Save?', "Are You Sure You Want To Save The Conditions?", save.Yes | save.No)
193         if user_click == save.Yes:
194             data={}
195             data['start_123'] = []
196             data['start_123'].append({
197                 'angle': int(self.start_123_angle.text()),
198                 'upper': int(self.start_123_upper.text()),
199                 'lower': int(self.start_123_lower.text())
200             })
201             data['start_234'] = []
202             data['start_234'].append({
203                 'angle': int(self.start_234_angle.text()),
204                 'upper': int(self.start_234_upper.text()),
205                 'lower': int(self.start_234_lower.text())
206             })
207
208             rep_data={}
209             rep_data['repetition']=[]
210             rep_data['repetition'].append({
211                 'repetition_number': int(self.Number_of_reps.value()),
212             })
213
214             with open(r"D:\STUDIES\YEAR 3 - SEM 1\ESP3902\openpose-master\build\examples\tutorial_api_python\Dev.Json\start.json", 'w') as outfile:
215                 json.dump(data, outfile)

```

```

216
217     with open(r"D:\STUDIES\YEAR 3 - SEM 1\ESP3902\openpose-master\build\examples\tutorial_api_python\Dev_Json\repetition.json", 'w') as outfile:
218         json.dump(rep_data, outfile)
219
220     saved = QMessageBox()
221     saved.setWindowTitle("Success")
222     saved.setText("Starting Pose Saved")
223     x = saved.exec_()
224
225
226 def save_end(self):
227     save = QMessageBox()
228     user_click = save.question(self,'End Pose Save?', "Are You Sure You Want To Save The Conditions?", save.Yes | save.No)
229     if(user_click == save.Yes):
230         data={}
231         data['end_234'] = []
232         data['end_234'].append({
233             'angle': int(self.end_234_angle.text()),
234             'upper': int(self.end_234_upper.text()),
235             'lower': int(self.end_234_lower.text())
236         })
237         data['end_814'] = []
238         data['end_814'].append({
239             'angle': int(self.end_814_angle.text()),
240             'upper': int(self.end_814_upper.text()),
241             'lower': int(self.end_814_lower.text())
242         })
243         data['end_13'] = []
244         data['end_13'].append({
245             'angle': int(self.end_13_angle.text()),
246             'upper': int(self.end_13_upper.text()),
247             'lower': int(self.end_13_lower.text())
248         })
249         data['end_14'] = []
250         data['end_14'].append({
251             'angle': int(self.end_14_angle.text()),
252             'upper': int(self.end_14_upper.text()),
253             'lower': int(self.end_14_lower.text())
254         })
255
256         rep_data={}
257         rep_data['repetition']=[]
258         rep_data['repetition'].append({
259             'repetition_number': int(self.Number_of_reps.value()),
260         })
261
262         with open(r"D:\STUDIES\YEAR 3 - SEM 1\ESP3902\openpose-master\build\examples\tutorial_api_python\Dev_Json\end.json", 'w') as outfile:
263             json.dump(data, outfile)
264
265         with open(r"D:\STUDIES\YEAR 3 - SEM 1\ESP3902\openpose-master\build\examples\tutorial_api_python\Dev_Json\repetition.json", 'w') as outfile:
266             json.dump(rep_data, outfile)
267
268         saved = QMessageBox()
269         saved.setWindowTitle("Success")
270         saved.setText("End Pose Saved")
271         x = saved.exec_()
272
273
274 def mousePressEvent(self, event):
275     self.dragPos = event.globalPos()
276
277
278
279
280
281
282 app = QtWidgets.QApplication(sys.argv)
283 window = Ui()
284 # app.exec_()
285 app.exec_()

```

Appendix II (Source code)

```
1 # Credit: https://github.com/intelRealSense/librealsense/issues/1904#issuecomment-398434434
2 import pyrealsense2 as rs
3 import numpy as np
4 import cv2
5 import os
6 import sys
7 from sys import platform
8 import argparse
9 import json
10 import math
11 import csv
12 import serial
13 from serial import *
14
15 # Realsense initialisation -----
16 # Configure depth and color streams
17 pipeline = rs.pipeline()
18
19 #Create a config and configure the pipeline to stream
20 config = rs.config()
21 config.enable_stream(rs.stream.depth, 640, 360, rs.format.z16, 30)
22 config.enable_stream(rs.stream.color, 640, 360, rs.format.bgr8, 30)
23
24 # Start streaming
25 pipe_profile = pipeline.start(config)
26
27 # Getting the depth sensor's depth scale (see rs-align example for explanation)
28 depth_sensor = pipe_profile.get_device().first_depth_sensor()
29 depth_scale = depth_sensor.get_depth_scale()
30
31 # Create an align object
32 # rs.align allows us to perform alignment of depth frames to others frames
33 # The "align_to" is the stream type to which we plan to align depth frames.
34 align_to = rs.stream.color
35 align = rs.align(align_to)
36 # RealSense end initialisation -----
37 |
38
39
40 # openpose initialisation -----
41 # Import Openpose
42 dir_path = os.path.dirname(os.path.realpath(__file__))
```

```

44 # Change these variables to point to the correct folder (Release/x64 etc.)
45 sys.path.append(dir_path + '../../../../../python/openpose/Release');
46 os.environ['PATH'] = os.environ['PATH'] + ';' + dir_path + '../../../../../x64/Release;' + dir_path + '../../../../../bin;';
47 import pyopenpose as op
48
49
50 # Flags
51 parser = argparse.ArgumentParser()
52 # parser.add_argument("--no_display", default=True, help="Enable to disable the visual display.")
53 args = parser.parse_known_args()
54
55 # !To use openpose in realsense, i need to be able to read image from a folder, then output the json to another folder,
56 # then realsense use that folder to import the json files
57 #
58 # Custom Params (refer to include/openpose/flags.hpp for more parameters)
59 params = dict()
60 params["model_folder"] = "../../../../../models/" # DO NOT CHANGE
61 params["display"] = 0
62 params["render_pose"] = 0
63 params["number_people_max"] = 1
64 params["write_json"] = "Test_input_images_json" # I think this is still needed. Can try removing it.
65
66
67 # Add others in path?
68 for i in range(0, len(args[1])):
69     curr_item = args[1][i]
70     if i != len(args[1])-1: next_item = args[1][i+1]
71     else: next_item = "1"
72     if "--" in curr_item and "--" in next_item:
73         key = curr_item.replace('--','')
74         if key not in params: params[key] = "1"
75     elif "--" in curr_item and "--" not in next_item:
76         key = curr_item.replace('--','')
77         if key not in params: params[key] = next_item
78
79 # Start openpose
80 opWrapper = opWrapperPython(0) # DONT TOUCH THE VALUE
81 opWrapper.configure(params)
82 opWrapper.start()
83 # openpose end initialisation -----
84
85
86
```

```

87     # Guidance skeleton initialisation -----
88     # END POSE
89     with open("D:/STUDIES/YEAR 3 - SEM 1/ESP3902/openpose-master/build/examples/tutorial_api_python/Pose Database/half_T_pose_end.json") as g:
90         guidance_skele_json = json.loads(g.read())
91
92     guidance_skele_part_end = guidance_skele_json["part_candidates"]
93
94     #STARTING POSE
95     with open("D:/STUDIES/YEAR 3 - SEM 1/ESP3902/openpose-master/build/examples/tutorial_api_python/Pose Database/half_T_pose_start.json") as g:
96         guidance_skele_json = json.loads(g.read())
97
98     guidance_skele_part_start = guidance_skele_json["part_candidates"]
99     # Guidance skeleton end initialisation -----
100
101
102
103     # MATH initialisation for GUIDANCE angles -----
104
105    with open(r"D:\STUDIES\YEAR 3 - SEM 1\ESP3902\openpose-master\build\examples\tutorial_api_python\Dev_Json\start.json") as start_json:
106        start = json.load(start_json)
107        for start_123 in start['start_123']:
108            guidance_start_angle123 = start_123["angle"]
109            start_angle123_upperAllowance = start_123["upper"]
110            start_angle123_lowerAllowance = start_123["lower"]
111        for start_234 in start['start_234']:
112            guidance_start_angle234 = start_234["angle"]
113            start_angle234_allowance = start_234["lower"]
114
115    with open(r"D:\STUDIES\YEAR 3 - SEM 1\ESP3902\openpose-master\build\examples\tutorial_api_python\Dev_Json\end.json") as end_json:
116        end = json.load(end_json)
117        for end_234 in end['end_234']:
118            guidance_end_angle234 = end_234["angle"]
119            end_angle234_allowance = end_234["lower"]
120        for end_814 in end['end_814']:
121            guidance_end_angle814 = end_814["angle"]
122            end_angle814_upperAllowance = end_814["upper"]
123            end_angle814_lowerAllowance = end_814["lower"]
124        for end_14 in end['end_14']:
125            guidance_planar_angle_chest_14 = end_14["angle"]
126            planar_angle_chest_14_upperAllowance = end_14["upper"]
127            planar_angle_chest_14_lowerAllowance = end_14["lower"]
128        for end_13 in end['end_13']:
129            guidance_planar_angle_chest_13 = end_13["angle"]

```

```

130     planar_angle_chest_13_upperAllowance = end_13["upper"]
131     planar_angle_chest_13_lowerAllowance = end_13["lower"]
132
133 with open(r"D:\STUDIES\YEAR 3 - SEM 1\ESP3902\openpose-master\build\examples\tutorial_api_python\Dev.Json\repetition.json") as rep_json:
134     rep = json.load(rep_json)
135     rep_value = rep['repetition'][0]['repetition_number']
136 # MATH end initialisation for GUIDANCE angles -----
137
138
139
140 # Arduino read from serial initialisation -----
141 ser = serial.Serial('COM7', 9600) # RMB TO CHANGE THE COM PORT BASE ON YOUR SETUP!
142 ser.xonoff=1
143 ser.flushInput()
144 # Arduino read from serial initialisation end -----
145
146
147
148 # Drawing config ----- *Realsense use BGR instead of RGB like a normal person. weirdo.
149 line_thickness = 2
150 colour_red = (0, 0, 255)
151 colour_green = (0, 128, 0)
152 colour_blue = (255,0,0)
153 font = cv2.FONT_HERSHEY_TRIPLEX
154 #-----
155
156
157 # global variables
158 curr_frame = 0
159 curr_pose = 0 # 0 is start, 1 is end
160 csv_frames = 0 # for CSV only
161 curr_rot = 0 # 0 is starting rotation, 1 is ending rotation part one, 2 is ending rotation part two
162 counter = 0 # for how many cycles completed
163 help_switch = 0 # the help screen is not toggled
164 bol_start = True # for the counter. Only after started then can increment.
165
166
167 # Streaming loop
168 try:
169     while True:
170         # Wait for a coherent pair of frames: depth and color
171         frames = pipeline.wait_for_frames() # frames.get_depth_frame() is a 640x360 depth image
172         depth_frame = frames.get_depth_frame()

```

```

173     color_frame = frames.get_color_frame()
174
175     # Align the depth frame to color frame
176     aligned_frames = align.process(frames)
177
178     # Get aligned frames
179     aligned_depth_frame = aligned_frames.get_depth_frame() # aligned_depth_frame is a 640x360 depth image
180     color_frame = aligned_frames.get_color_frame()
181
182     # Validate that both frames are valid
183     if not aligned_depth_frame or not color_frame:
184         continue
185
186     # Intrinsic & Extrinsic
187     depth_intrinsic = aligned_depth_frame.profile.as_video_stream_profile().intrinsic
188     color_intrinsic = color_frame.profile.as_video_stream_profile().intrinsic
189     depth_to_color_extrinsic = aligned_depth_frame.profile.get_extrinsics_to(color_frame.profile)
190
191     # Convert images to numpy arrays
192     depth_image = np.asarray(aligned_depth_frame.get_data())
193     color_image = np.asarray(color_frame.get_data())
194
195     # save the colour image into a folder in path
196     path = "D:/STUDIES/YEAR 3 - SEM 1/ESP3902/openpose-master/build/examples/tutorial_api_python/Realsense_colour_images"
197     cv2.imwrite(os.path.join(path , str(curr_frame) + '_colour_frame.png'), color_image)
198
199     # OpenPose
200     imagepath = r'D:/STUDIES/YEAR 3 - SEM 1/ESP3902/openpose-master/build/examples/tutorial_api_python/Realsense_colour_images\0.colour_frame.png'
201     datum = op.Datum()
202     imageToProcess = cv2.imread(imagepath)
203     datum.cvInputData = imageToProcess
204     opWrapper.emplaceAndPop([datum])
205     np.userskele_json_part = datum.poseKeypoints
206     # if(np.userskele_json_part.ndim > 0): # !!IMPORTANT!! MAKE SURE THERE IS THIS CONDITION TO MAKE THE VARIABLE HAVE PROPER KEYPOINT
207     #     print(np.userskele_json_part[0][0][0])
208     # if(datum.poseKeypoints.ndim > 0): # not having to copy the array over to a new variable works as well
209     #     print(datum.poseKeypoints[0][4])
210
211
212     # Stack both images horizontally, but now i only want the colour_image
213     images = color_image
214
215 #####

```

```

216
217     # user skeleton calculation # NEW VERSION (DATUM)
218
219     if(datum.poseKeypoints.ndim > 0):
220         # 123
221         if(datum.poseKeypoints[0][2].all and datum.poseKeypoints[0][1].all and datum.poseKeypoints[0][3].all):
222             if(any(datum.poseKeypoints[0][2]) and any(datum.poseKeypoints[0][1]) and any(datum.poseKeypoints[0][3])):
223                 # print("\n123")
224                 AB = [datum.poseKeypoints[0][1][0] - datum.poseKeypoints[0][2][0], datum.poseKeypoints[0][1][1] - datum.poseKeypoints[0][2][1]]
225                 AC = [datum.poseKeypoints[0][3][0] - datum.poseKeypoints[0][2][0], datum.poseKeypoints[0][3][1] - datum.poseKeypoints[0][2][1]]
226                 # dot_prdt = np.array(AB) @ np.array(AC)
227                 dot_prdt = np.dot(AB, AC)
228                 # print("dot_prdt: ", dot_prdt)
229                 len_AB = math.sqrt(AB[0]**2 + AB[1]**2)
230                 len_AC = math.sqrt(AC[0]**2 + AC[1]**2)
231                 Length = len_AB * len_AC
232                 # print("Length: ", Length)
233                 if(len_AB!=0 and len_AC!=0):
234                     divided_value = np.divide(dot_prdt,Length)
235                     # print("divided_value: ", divided_value)
236                     if(divided_value<1 and divided_value>-1):
237                         user_angle123 = math.degrees(math.acos(divided_value))
238                         # print("user_angle123: ", user_angle123)
239                     else:
240                         user_angle123=0
241                 else:
242                     user_angle123 = 0
243             user_angle123 = 0 # needed else will not work. So in case when either 1 of the parts is null, then just let the angle =0
244
245             # 234
246             if(datum.poseKeypoints[0][2].all and datum.poseKeypoints[0][3].all and datum.poseKeypoints[0][4].all):
247                 if(any(datum.poseKeypoints[0][2]) and any(datum.poseKeypoints[0][3]) and any(datum.poseKeypoints[0][4])):
248                     # print("\n234")
249                     AB = [datum.poseKeypoints[0][2][0] - datum.poseKeypoints[0][3][0], datum.poseKeypoints[0][2][1] - datum.poseKeypoints[0][3][1]]
250                     AC = [datum.poseKeypoints[0][4][0] - datum.poseKeypoints[0][3][0], datum.poseKeypoints[0][4][1] - datum.poseKeypoints[0][3][1]]
251                     # dot_prdt = np.array(AB) @ np.array(AC)
252                     dot_prdt = np.dot(AB, AC)
253                     # print("dot_prdt: ", dot_prdt)
254                     len_AB = math.sqrt(AB[0]**2 + AB[1]**2)
255                     len_AC = math.sqrt(AC[0]**2 + AC[1]**2)
256                     Length = len_AB * len_AC
257                     # print("Length: ", Length)
258                     if(len_AB!=0 and len_AC!=0):

```

```

259         divided_value = np.divide(dot_prdt,length)
260         # print("divided_value: ", divided_value)
261         if(divided_value<1 and divided_value>-1):
262             user_angle234 = math.degrees(math.acos(divided_value))
263             user_angle234 = math.degrees(math.acos(divided_value))
264             # print("user_angle234: ", user_angle234)
265         else:
266             user_angle234=0
267         else:
268             user_angle234=0
269
270         # print(user_angle234)
271     else:
272         user_angle234 = 0
273
274     user_angle234 = 0 # needed else will not work. So in case when either 1 of the parts is null, then just let the angle =0
275
276 # 814
277 if(datum.poseKeypoints[0][8].all and datum.poseKeypoints[0][1].all and datum.poseKeypoints[0][4].all):
278     if(any(datum.poseKeypoints[0][2]) and any(datum.poseKeypoints[0][1]) and any(datum.poseKeypoints[0][4])):
279         # print("\n813")
280         AB = [datum.poseKeypoints[0][8][0] - datum.poseKeypoints[0][1][0], datum.poseKeypoints[0][8][1] - datum.poseKeypoints[0][1][1]]
281         AC = [datum.poseKeypoints[0][4][0] - datum.poseKeypoints[0][1][0], datum.poseKeypoints[0][4][1] - datum.poseKeypoints[0][1][1]]
282         # dot_prdt = np.array(AB) @ np.array(AC)
283         dot_prdt = np.dot(AB, AC)
284         # print("dot_prdt: ", dot_prdt)
285         len_AB = math.sqrt(AB[0]**2 + AB[1]**2)
286         len_AC = math.sqrt(AC[0]**2 + AC[1]**2)
287         Length = len_AB * len_AC
288         # print("Length: ", Length)
289         if(len_AB!=0 and len_AC!=0):
290             divided_value = np.divide(dot_prdt,Length)
291             # print("divided_value: ", divided_value)
292             if(divided_value<1 and divided_value>-1):
293                 user_angle814 = math.degrees(math.acos(divided_value))
294                 # print("user_angle814: ", user_angle814)
295             else:
296                 user_angle814=0
297             else:
298                 user_angle814 = 0
299
300     user_angle814 = 0 # needed else will not work. So in case when either 1 of the parts is null, then just let the angle =0
301
302 # Normal of 1-9-12 @ 1-4 @ 1-3
303 if(datum.poseKeypoints[0][1].all and datum.poseKeypoints[0][9].all and datum.poseKeypoints[0][12].all and datum.poseKeypoints[0][4].all and datum.poseKeypoints[0][3].all):
304     if((any(datum.poseKeypoints[0][1]) and any(datum.poseKeypoints[0][9])) and any(datum.poseKeypoints[0][12]) and any(datum.poseKeypoints[0][4]) and any(datum.poseKeypoints[0][3])):
305         depth_keypoint1 = aligned_depth_frame.get_distance(int(datum.poseKeypoints[0][1][0]), int(datum.poseKeypoints[0][1][1])) #A
306         depth_keypoint9 = aligned_depth_frame.get_distance(int(datum.poseKeypoints[0][9][0]), int(datum.poseKeypoints[0][9][1])) #B
307         depth_keypoint12 = aligned_depth_frame.get_distance(int(datum.poseKeypoints[0][12][0]), int(datum.poseKeypoints[0][12][1])) #C
308         depth_keypoint3 = aligned_depth_frame.get_distance(int(datum.poseKeypoints[0][3][0]), int(datum.poseKeypoints[0][3][1])) #D
309         depth_keypoint4 = aligned_depth_frame.get_distance(int(datum.poseKeypoints[0][4][0]), int(datum.poseKeypoints[0][4][1])) #E
310
311         point1 = rs.rs2_deproject_pixel_to_point(color_intrin, [int(datum.poseKeypoints[0][1][0]), int(datum.poseKeypoints[0][1][1])], depth_keypoint1)
312         point9 = rs.rs2_deproject_pixel_to_point(color_intrin, [int(datum.poseKeypoints[0][9][0]), int(datum.poseKeypoints[0][9][1])], depth_keypoint9)
313         point12 = rs.rs2_deproject_pixel_to_point(color_intrin, [int(datum.poseKeypoints[0][12][0]), int(datum.poseKeypoints[0][12][1])], depth_keypoint12)
314         point4 = rs.rs2_deproject_pixel_to_point(color_intrin, [int(datum.poseKeypoints[0][4][0]), int(datum.poseKeypoints[0][4][1])], depth_keypoint4)
315         point3 = rs.rs2_deproject_pixel_to_point(color_intrin, [int(datum.poseKeypoints[0][3][0]), int(datum.poseKeypoints[0][3][1])], depth_keypoint3)
316
317         len_AB = math.sqrt(math.pow(point1[0] - point9[0], 2) + math.pow(point1[1] - point9[1], 2) + math.pow(point1[2] - point9[2], 2)) #1-9
318         len_AC = math.sqrt(math.pow(point1[0] - point12[0], 2) + math.pow(point1[1] - point12[1], 2) + math.pow(point1[2] - point12[2], 2)) #1-12
319         len_AE = math.sqrt(math.pow(point1[0] - point4[0], 2) + math.pow(point1[1] - point4[1], 2) + math.pow(point1[2] - point4[2], 2)) #1-4
320         len_AD = math.sqrt(math.pow(point1[0] - point3[0], 2) + math.pow(point1[1] - point3[1], 2) + math.pow(point1[2] - point3[2], 2)) #1-3
321
322         AB = [point9[0] - point1[0], point9[1] - point1[1], depth_keypoint9-depth_keypoint1]
323         AC = [point12[0] - point1[0], point12[1] - point1[1], depth_keypoint12-depth_keypoint1]
324         AE = [point4[0] - point1[0], point4[1] - point1[1], depth_keypoint4-depth_keypoint1]
325         AD = [point3[0] - point1[0], point3[1] - point1[1], depth_keypoint3-depth_keypoint1]
326
327         if(len_AB != 0 and len_AC != 0 and len_AE != 0 and len_AD != 0):
328             cross_prdt = np.cross(AB, AC)
329             len_cross_prdt=math.sqrt(math.pow(cross_prdt[0], 2) + math.pow(cross_prdt[1], 2) + math.pow(cross_prdt[2], 2))
330             norm_cross_prdt = np.divide(cross_prdt, len_cross_prdt)
331             norm_dot_prdt = np.divide((np.array(norm_cross_prdt) @ np.array(AE)), (len_AE))
332             planar_angle_chest_14 = math.degrees(math.acos(norm_dot_prdt))
333             norm_dot_prdt = np.divide((np.array(norm_cross_prdt) @ np.array(AD)), (len_AD))
334             planar_angle_chest_13 = math.degrees(math.acos(norm_dot_prdt))
335
336             # print("planar: ", planar_angle_chest)
337
338             # print(planar_angle_chest)
339
340             # distance away from camera warning
341             # print(depth_keypoint1)
342             if(depth_keypoint1>3.0):
343                 images = cv2.rectangle(images, (5, 340),(140, 360), (255, 255, 255), -1)
344                 cv2.putText(images,'Move forward',(7, 355),font, 0.5,(0, 255, 255),cv2.LINE_4)

```

```

345
346         bol_wrong = True
347     elif(depth_keypoint1<2.0):
348         # print("AJKLSDHLAJSDLKJASLKDJALKSJDLKAJSSDLKJASLKDJALKSJDLKAJSSDLKJASLDKJASe")
349         images = cv2.rectangle(images, (5, 340) ,(140, 360) , (255, 255, 255), -1)
350         cv2.putText(images,"Move backward", (7, 355),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
351         # print("asdasdasdasdasdasdasdasdas")
352         bol_wrong = True
353     else:
354         bol_wrong = False
355
356     user_angle123 = 0
357     user_angle234 = 0
358     user_angle814 = 0
359     planar_angle_chest_14 = 0
360     planar_angle_chest_13 = 0
361
362     # print("123:",user_angle123)
363     # print("234:",user_angle234)
364     # print("814:",user_angle814)
365     # print("14:",planar_angle_chest_14)
366     # print("13:",planar_angle_chest_13)
367
368     # User rotation axis. 0-2 in terms of x,y,z
369     ser_bytes = ser.readline()
370     # print(ser_bytes)
371     decoded_bytes = ser_bytes[0:len(ser_bytes)].decode("utf-8").strip('\n').split(',')
372     # print(decoded_bytes)
373     if"\x00" not in decoded_bytes[0] and decoded_bytes[0]!='':
374         # if decoded_bytes[0].isdecimal() == True:
375             float_list_values = [float(i) for i in decoded_bytes]
376             highest_list_values = max(float_list_values, key=abs)
377             highest_axis_index = float_list_values.index(highest_list_values)
378             # print(highest_list_values)
379             # print("rot index:",highest_axis_index)
380     else:
381         highest_axis_index = 2 # since we not using z-axis
382         highest_list_values = 0 # doesnt matter what i put
383     # problem with not reading the latest data from arduino. it reads the top in a ever-filling container of data from arduino.
384     # tried delay => fps worst
385     # tried https://stackoverflow.com/questions/1093598/pyserial-how-to-read-the-last-line-sent-from-a-serial-device => doesnt work
386     # tried https://stackoverflow.com/questions/1093598/pyserial-how-to-read-the-last-line-sent-from-a-serial-device => doesnt work
387
388
389 # User skeleton drawing # NEW VERSION (DATUM)
390 if(datum.posekeypoints.ndim > 0):
391     for i in range(15):
392         if(datum.posekeypoints[0][i].all): # need to .all to make sure that that keypoint returns true for it to be drawn
393             if(any(datum.posekeypoints[0][i])):
394                 cv2.circle(images, (int(datum.posekeypoints[0][i][0]),int(datum.posekeypoints[0][i][1])), 5, colour_blue, thickness=line_thickness)
395
396
397 # Guidance skeleton drawing (Line version)
398 if(curr_pose == 1): # change to end pose
399     if((user_angle234>guidance_end_angle234_end_angle234_allowance) or (user_angle814>guidance_end_angle814_end_angle814_lowerAllowance or user_angle814>guidance_end_angle814_end_angle814_upperAllowance) or
400         (planar_angle_chest_13>guidance_planar_angle_chest_13_planar_angle_chest_13_lowerAllowance or planar_angle_chest_13>guidance_planar_angle_chest_13_planar_angle_chest_13_upperAllowance) or
401         (planar_angle_chest_14>guidance_planar_angle_chest_14_planar_angle_chest_14_lowerAllowance or planar_angle_chest_14>guidance_planar_angle_chest_14_planar_angle_chest_14_upperAllowance) or
402         (curr_rot==1 and (highest_axis_index!=0 or highest_list_values<0)) or (curr_rot==2 and (highest_axis_index!=0 or highest_list_values>0)) or
403         bol_wrong == True): # outside tolerance, red colour
404         # 0 - 1
405         cv2.line(images, (int(guidance_skele_part_end[0]['0'][0]),int(guidance_skele_part_end[0]['0'][1])), (int(guidance_skele_part_end[0]['1'][0]),int(guidance_skele_part_end[0]['1'][1])), colour_red, thickness=line_thickness)
406         # 1 - 2
407         cv2.line(images, (int(guidance_skele_part_end[0]['2'][0]),int(guidance_skele_part_end[0]['2'][1])), (int(guidance_skele_part_end[0]['1'][0]),int(guidance_skele_part_end[0]['1'][1])), colour_red, thickness=line_thickness)
408         # 2 - 3
409         cv2.line(images, (int(guidance_skele_part_end[0]['5'][0]),int(guidance_skele_part_end[0]['5'][1])), (int(guidance_skele_part_end[0]['1'][0]),int(guidance_skele_part_end[0]['1'][1])), colour_red, thickness=line_thickness)
410         # 3 - 8
411         cv2.line(images, (int(guidance_skele_part_end[0]['8'][0]),int(guidance_skele_part_end[0]['8'][1])), (int(guidance_skele_part_end[0]['1'][0]),int(guidance_skele_part_end[0]['1'][1])), colour_red, thickness=line_thickness)
412         # 2 - 3
413         cv2.line(images, (int(guidance_skele_part_end[0]['2'][0]),int(guidance_skele_part_end[0]['2'][1])), (int(guidance_skele_part_end[0]['3'][0]),int(guidance_skele_part_end[0]['3'][1])), colour_red, thickness=line_thickness)
414         # 3 - 4
415         cv2.line(images, (int(guidance_skele_part_end[0]['3'][0]),int(guidance_skele_part_end[0]['3'][1])), (int(guidance_skele_part_end[0]['4'][0]),int(guidance_skele_part_end[0]['4'][1])), colour_red, thickness=line_thickness)
416         # 5 - guidance_skele_part_end
417         cv2.line(images, (int(guidance_skele_part_end[0]['5'][0]),int(guidance_skele_part_end[0]['5'][1])), (int(guidance_skele_part_end[0]['6'][0]),int(guidance_skele_part_end[0]['6'][1])), colour_red, thickness=line_thickness)
418         # 6 - 7
419         cv2.line(images, (int(guidance_skele_part_end[0]['6'][0]),int(guidance_skele_part_end[0]['6'][1])), (int(guidance_skele_part_end[0]['7'][0]),int(guidance_skele_part_end[0]['7'][1])), colour_red, thickness=line_thickness)
420         # 8 - 9
421         cv2.line(images, (int(guidance_skele_part_end[0]['8'][0]),int(guidance_skele_part_end[0]['8'][1])), (int(guidance_skele_part_end[0]['9'][0]),int(guidance_skele_part_end[0]['9'][1])), colour_red, thickness=line_thickness)
422         # 8 - 12
423         cv2.line(images, (int(guidance_skele_part_end[0]['8'][0]),int(guidance_skele_part_end[0]['8'][1])), (int(guidance_skele_part_end[0]['12'][0]),int(guidance_skele_part_end[0]['12'][1])), colour_red, thickness=line_thickness)
424         # 9 - 10
425         cv2.line(images, (int(guidance_skele_part_end[0]['9'][0]),int(guidance_skele_part_end[0]['9'][1])), (int(guidance_skele_part_end[0]['10'][0]),int(guidance_skele_part_end[0]['10'][1])), colour_red, thickness=line_thickness)
426         # 10 - 11
427         cv2.line(images, (int(guidance_skele_part_end[0]['10'][0]),int(guidance_skele_part_end[0]['10'][1])), (int(guidance_skele_part_end[0]['11'][0]),int(guidance_skele_part_end[0]['11'][1])), colour_red, thickness=line_thickness)
428         # 12 - 13
429         cv2.line(images, (int(guidance_skele_part_end[0]['12'][0]),int(guidance_skele_part_end[0]['12'][1])), (int(guidance_skele_part_end[0]['13'][0]),int(guidance_skele_part_end[0]['13'][1])), colour_red, thickness=line_thickness)
430         # 13 - 14
431         cv2.line(images, (int(guidance_skele_part_end[0]['13'][0]),int(guidance_skele_part_end[0]['13'][1])), (int(guidance_skele_part_end[0]['14'][0]),int(guidance_skele_part_end[0]['14'][1])), colour_red, thickness=line_thickness)
432
433     # Error warnings
434     #rotation
435     if(curr_rot==1 and (highest_axis_index!=0 or highest_list_values<0)):
436         # Draw a rectangle
437         images = cv2.rectangle(images, (0, 0), (213, 70), (255, 255, 255), -1) # frame to draw, top left, bot right, colour, -1
438         cv2.putText(images,'Rotate your hand until',(5, 15),font, 0.5,(0, 0, 255),1,cv2.LINE_4)

```

```

439         cv2.putText(images,'your palm faces the',(5, 40),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
440         cv2.putText(images,'floor',(5, 65),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
441     elif(curr_rot==2 and (highest_axis_index!=0 or highest_list_values!=0)):
442         # Draw a rectangle
443         images = cv2.rectangle(images, (0, 0), (213, 70), (255, 255, 255), -1)
444         cv2.putText(images,'Rotate your hand until',(5, 15),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
445         cv2.putText(images,'your palm faces the',(5, 40),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
446         cv2.putText(images,'ceiling',(5, 65),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
447     # elbow
448     if(user_angle234==guidance_end_angle234-end_angle234_allowance):
449         images = cv2.rectangle(images, (225, 340) , (425, 360) , (255, 255, 255), -1)
450         cv2.putText(images,'Straighten your arm',(230, 355),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
451     # up, down
452     if(user_angle814==guidance_end_angle814-end_angle814_lowallowance):
453         images = cv2.rectangle(images, (228, 0) , (405, 20) , (255, 255, 255), -1)
454         cv2.putText(images,'Move your arm up',(232, 15),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
455     elif(user_angle814==guidance_end_angle814-end_angle814_upperallowance):
456         images = cv2.rectangle(images, (216, 0) , (406, 20) , (255, 255, 255), -1)
457         cv2.putText(images,'Move your arm down',(218, 15),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
458     # front, back
459     if(planar_angle_chest_13==guidance_planar_angle_chest_13_lowallowance and planar_angle_chest_14==guidance_planar_angle_chest_14_lowallowance):
460         images = cv2.rectangle(images, (410, 0) ,(640, 20) , (255, 255, 255), -1)
461         cv2.putText(images,'Move your arm backward',(413, 15),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
462     elif(planar_angle_chest_13==guidance_planar_angle_chest_13_upperallowance and planar_angle_chest_14==guidance_planar_angle_chest_14_upperallowance):
463         images = cv2.rectangle(images, (420, 0) ,(640, 20) , (255, 255, 255), -1)
464         cv2.putText(images,'Move your arm forward',(425, 15),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
465     else:
466         images = cv2.rectangle(images, (420, 0) ,(640, 20) , (255, 255, 255), -1)
467         cv2.putText(images,'Straighten your elbow',(425, 15),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
468
469 else: # within tolerance, green colour
470     ## - 1
471     cv2.line(images, (int(guidance_skele_part_end[0]['0'][0]),int(guidance_skele_part_end[0]['0'][1])), (int(guidance_skele_part_end[0]['1'][0]),int(guidance_skele_part_end[0]['1'][1])), colour_green, thickness=line_thickness)
472     ## - 2
473     cv2.line(images, (int(guidance_skele_part_end[0]['2'][0]),int(guidance_skele_part_end[0]['2'][1])), (int(guidance_skele_part_end[0]['1'][0]),int(guidance_skele_part_end[0]['1'][1])), colour_green, thickness=line_thickness)
474     ## - 5
475     cv2.line(images, (int(guidance_skele_part_end[0]['5'][0]),int(guidance_skele_part_end[0]['5'][1])), (int(guidance_skele_part_end[0]['1'][0]),int(guidance_skele_part_end[0]['1'][1])), colour_green, thickness=line_thickness)
476     ## - 8
477     cv2.line(images, (int(guidance_skele_part_end[0]['8'][0]),int(guidance_skele_part_end[0]['8'][1])), (int(guidance_skele_part_end[0]['1'][0]),int(guidance_skele_part_end[0]['1'][1])), colour_green, thickness=line_thickness)
478     ## - 3
479     cv2.line(images, (int(guidance_skele_part_end[0]['2'][0]),int(guidance_skele_part_end[0]['2'][1])), (int(guidance_skele_part_end[0]['3'][0]),int(guidance_skele_part_end[0]['3'][1])), colour_green, thickness=line_thickness)
480     ## - 4
481     cv2.line(images, (int(guidance_skele_part_end[0]['3'][0]),int(guidance_skele_part_end[0]['3'][1])), (int(guidance_skele_part_end[0]['4'][0]),int(guidance_skele_part_end[0]['4'][1])), colour_green, thickness=line_thickness)
482     ## - guidance_skele_part_end
483     cv2.line(images, (int(guidance_skele_part_end[0]['5'][0]),int(guidance_skele_part_end[0]['5'][1])), (int(guidance_skele_part_end[0]['6'][0]),int(guidance_skele_part_end[0]['6'][1])), colour_green, thickness=line_thickness)
484     ## - 7
485     cv2.line(images, (int(guidance_skele_part_end[0]['6'][0]),int(guidance_skele_part_end[0]['6'][1])), (int(guidance_skele_part_end[0]['7'][0]),int(guidance_skele_part_end[0]['7'][1])), colour_green, thickness=line_thickness)
486     ## - 9
487     cv2.line(images, (int(guidance_skele_part_end[0]['8'][0]),int(guidance_skele_part_end[0]['8'][1])), (int(guidance_skele_part_end[0]['9'][0]),int(guidance_skele_part_end[0]['9'][1])), colour_green, thickness=line_thickness)
488     ## - 12
489     cv2.line(images, (int(guidance_skele_part_end[0]['8'][0]),int(guidance_skele_part_end[0]['8'][1])), (int(guidance_skele_part_end[0]['12'][0]),int(guidance_skele_part_end[0]['12'][1])), colour_green, thickness=line_thickness)

490     cv2.line(images, (int(guidance_skele_part_end[0]['0'][0]),int(guidance_skele_part_end[0]['1'][0])), (int(guidance_skele_part_end[0]['12'][0]),int(guidance_skele_part_end[0]['12'][1])), colour_green, thickness=line_thickness)
491     ## - 10
492     cv2.line(images, (int(guidance_skele_part_end[0]['9'][0]),int(guidance_skele_part_end[0]['9'][1])), (int(guidance_skele_part_end[0]['10'][0]),int(guidance_skele_part_end[0]['10'][1])), colour_green, thickness=line_thickness)
493     ## - 11
494     cv2.line(images, (int(guidance_skele_part_end[0]['10'][0]),int(guidance_skele_part_end[0]['10'][1])), (int(guidance_skele_part_end[0]['11'][0]),int(guidance_skele_part_end[0]['11'][1])), colour_green, thickness=line_thickness)
495     ## - 13
496     cv2.line(images, (int(guidance_skele_part_end[0]['12'][0]),int(guidance_skele_part_end[0]['12'][1])), (int(guidance_skele_part_end[0]['13'][0]),int(guidance_skele_part_end[0]['13'][1])), colour_green, thickness=line_thickness)
497     ## - 14
498     cv2.line(images, (int(guidance_skele_part_end[0]['13'][0]),int(guidance_skele_part_end[0]['13'][1])), (int(guidance_skele_part_end[0]['14'][0]),int(guidance_skele_part_end[0]['14'][1])), colour_green, thickness=line_thickness)
499     #RENAME TO START POSE REFLICK PART TWO OF ROTATION
500     if(curr_ang == 1):
501         curr_rot=2
502     elif(curr_ang == 2):
503         curr_pose=0
504         curr_rot=0
505
506 else: # starting pose
507     if((user_angle123==guidance_start_angle123-start_angle123_lowallowance or user_angle123==guidance_start_angle123-start_angle123_upperallowance) or user_angle234==guidance_start_angle234-start_angle234_allowance) or (highest_axis_index!=1 or highest_list_values!=0) or
508     hol_wrong == True): # outside tolerance, red colour
509     ## - 1
510     cv2.line(images, (int(guidance_skele_part_start[0]['0'][0]),int(guidance_skele_part_start[0]['0'][1])), (int(guidance_skele_part_start[0]['1'][0]),int(guidance_skele_part_start[0]['1'][1])), colour_red, thickness=line_thickness)
511     ## - 2
512     cv2.line(images, (int(guidance_skele_part_start[0]['2'][0]),int(guidance_skele_part_start[0]['2'][1])), (int(guidance_skele_part_start[0]['1'][0]),int(guidance_skele_part_start[0]['1'][1])), colour_red, thickness=line_thickness)
513     ## - 5
514     cv2.line(images, (int(guidance_skele_part_start[0]['5'][0]),int(guidance_skele_part_start[0]['5'][1])), (int(guidance_skele_part_start[0]['1'][0]),int(guidance_skele_part_start[0]['1'][1])), colour_red, thickness=line_thickness)
515     ## - 8
516     cv2.line(images, (int(guidance_skele_part_start[0]['8'][0]),int(guidance_skele_part_start[0]['8'][1])), (int(guidance_skele_part_start[0]['1'][0]),int(guidance_skele_part_start[0]['1'][1])), colour_red, thickness=line_thickness)
517     ## - 3
518     cv2.line(images, (int(guidance_skele_part_start[0]['2'][0]),int(guidance_skele_part_start[0]['2'][1])), (int(guidance_skele_part_start[0]['3'][0]),int(guidance_skele_part_start[0]['3'][1])), colour_red, thickness=line_thickness)
519     ## - 4
520     cv2.line(images, (int(guidance_skele_part_start[0]['3'][0]),int(guidance_skele_part_start[0]['3'][1])), (int(guidance_skele_part_start[0]['4'][0]),int(guidance_skele_part_start[0]['4'][1])), colour_red, thickness=line_thickness)
521     ## - 6
522     cv2.line(images, (int(guidance_skele_part_start[0]['5'][0]),int(guidance_skele_part_start[0]['5'][1])), (int(guidance_skele_part_start[0]['6'][0]),int(guidance_skele_part_start[0]['6'][1])), colour_red, thickness=line_thickness)
523     ## - 8
524     cv2.line(images, (int(guidance_skele_part_start[0]['6'][0]),int(guidance_skele_part_start[0]['6'][1])), (int(guidance_skele_part_start[0]['7'][0]),int(guidance_skele_part_start[0]['7'][1])), colour_red, thickness=line_thickness)
525     ## - 13
526     cv2.line(images, (int(guidance_skele_part_start[0]['12'][0]),int(guidance_skele_part_start[0]['12'][1])), (int(guidance_skele_part_start[0]['13'][0]),int(guidance_skele_part_start[0]['13'][1])), colour_red, thickness=line_thickness)
527     ## - 12
528     cv2.line(images, (int(guidance_skele_part_start[0]['14'][0]),int(guidance_skele_part_start[0]['14'][1])), (int(guidance_skele_part_start[0]['12'][0]),int(guidance_skele_part_start[0]['12'][1])), colour_red, thickness=line_thickness)
529     ## - 10
530     cv2.line(images, (int(guidance_skele_part_start[0]['9'][0]),int(guidance_skele_part_start[0]['9'][1])), (int(guidance_skele_part_start[0]['10'][0]),int(guidance_skele_part_start[0]['10'][1])), colour_red, thickness=line_thickness)
531     ## - 11
532     cv2.line(images, (int(guidance_skele_part_start[0]['10'][0]),int(guidance_skele_part_start[0]['10'][1])), (int(guidance_skele_part_start[0]['11'][0]),int(guidance_skele_part_start[0]['11'][1])), colour_red, thickness=line_thickness)
533     ## - 13
534     cv2.line(images, (int(guidance_skele_part_start[0]['12'][0]),int(guidance_skele_part_start[0]['12'][1])), (int(guidance_skele_part_start[0]['13'][0]),int(guidance_skele_part_start[0]['13'][1])), colour_red, thickness=line_thickness)
535     ## - 14
536     cv2.line(images, (int(guidance_skele_part_start[0]['13'][0]),int(guidance_skele_part_start[0]['13'][1])), (int(guidance_skele_part_start[0]['14'][0]),int(guidance_skele_part_start[0]['14'][1])), colour_red, thickness=line_thickness)
537
538     #rotation
539     if(highest_axis_index!=1 or highest_list_values!=0):
540         images = cv2.rectangle(images, (0, 0) ,(235, 50), (255, 255, 255), -1) # Frame to draw, top left, bot right, colour, -1
541         cv2.putText(images,'Rotate your hand until',(5, 15),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
542         cv2.putText(images,'your palm faces your leg',(5, 40),font, 0.5,(0, 0, 255),1,cv2.LINE_4)
543
544     else: # within tolerance, green colour
545     ## - 1
546     cv2.line(images, (int(guidance_skele_part_start[0]['0'][0]),int(guidance_skele_part_start[0]['0'][1])), (int(guidance_skele_part_start[0]['1'][0]),int(guidance_skele_part_start[0]['1'][1])), colour_green, thickness=line_thickness)
547     ## - 2
548     cv2.line(images, (int(guidance_skele_part_start[0]['2'][0]),int(guidance_skele_part_start[0]['2'][1])), (int(guidance_skele_part_start[0]['1'][0]),int(guidance_skele_part_start[0]['1'][1])), colour_green, thickness=line_thickness)

```

```

551    #1 - 5
552    cv2.line(images, (int(guidance_skele_part_start[0]['5'][0]),int(guidance_skele_part_start[0]['5'][1])), (int(guidance_skele_part_start[0]['1'][0]),int(guidance_skele_part_start[0]['1'][1])), colour_green, thickness=line_thickness)
553    #1 - 6
554    cv2.line(images, (int(guidance_skele_part_start[0]['6'][0]),int(guidance_skele_part_start[0]['6'][1])), (int(guidance_skele_part_start[0]['1'][0]),int(guidance_skele_part_start[0]['1'][1])), colour_green, thickness=line_thickness)
555    #2 - 3
556    cv2.line(images, (int(guidance_skele_part_start[0]['2'][0]),int(guidance_skele_part_start[0]['2'][1])), (int(guidance_skele_part_start[0]['3'][0]),int(guidance_skele_part_start[0]['3'][1])), colour_green, thickness=line_thickness)
557    #3 - 4
558    cv2.line(images, (int(guidance_skele_part_start[0]['4'][0]),int(guidance_skele_part_start[0]['3'][1])), (int(guidance_skele_part_start[0]['4'][0]),int(guidance_skele_part_start[0]['3'][1])), colour_green, thickness=line_thickness)
559    #5 - 6
560    cv2.line(images, (int(guidance_skele_part_start[0]['6'][0]),int(guidance_skele_part_start[0]['5'][1])), (int(guidance_skele_part_start[0]['6'][0]),int(guidance_skele_part_start[0]['6'][1])), colour_green, thickness=line_thickness)
561    #6 - 7
562    cv2.line(images, (int(guidance_skele_part_start[0]['7'][0]),int(guidance_skele_part_start[0]['6'][1])), (int(guidance_skele_part_start[0]['7'][0]),int(guidance_skele_part_start[0]['7'][1])), colour_green, thickness=line_thickness)
563    #8 - 9
564    cv2.line(images, (int(guidance_skele_part_start[0]['9'][0]),int(guidance_skele_part_start[0]['8'][1])), (int(guidance_skele_part_start[0]['9'][0]),int(guidance_skele_part_start[0]['9'][1])), colour_green, thickness=line_thickness)
565    #10 - 12
566    cv2.line(images, (int(guidance_skele_part_start[0]['12'][0]),int(guidance_skele_part_start[0]['10'][1])), (int(guidance_skele_part_start[0]['12'][0]),int(guidance_skele_part_start[0]['12'][1])), colour_green, thickness=line_thickness)
567    #9 - 10
568    cv2.line(images, (int(guidance_skele_part_start[0]['10'][0]),int(guidance_skele_part_start[0]['9'][1])), (int(guidance_skele_part_start[0]['10'][0]),int(guidance_skele_part_start[0]['10'][1])), colour_green, thickness=line_thickness)
569    #11 - 12
570    cv2.line(images, (int(guidance_skele_part_start[0]['10'][0]),int(guidance_skele_part_start[0]['11'][1])), (int(guidance_skele_part_start[0]['11'][0]),int(guidance_skele_part_start[0]['11'][1])), colour_green, thickness=line_thickness)
571    #12 - 13
572    cv2.line(images, (int(guidance_skele_part_start[0]['12'][0]),int(guidance_skele_part_start[0]['12'][1])), (int(guidance_skele_part_start[0]['13'][0]),int(guidance_skele_part_start[0]['13'][1])), colour_green, thickness=line_thickness)
573    #13 - 14
574    cv2.line(images, (int(guidance_skele_part_start[0]['13'][0]),int(guidance_skele_part_start[0]['13'][1])), (int(guidance_skele_part_start[0]['14'][0]),int(guidance_skele_part_start[0]['14'][1])), colour_green, thickness=line_thickness)
575    #CHANGE TO END POSE
576    CURT_pose = 1
577    curt_rot = 1
578    if bol_start != True:
579        counter += 1
580    bol_start = False
581
582 csv_frames += 1
583
584 # print("User_angle123a ", user_angle123a)
585 # print("User_angle1234 ", user_angle1234)
586 # print("User_angle101a ", user_angle101a)
587 # print("planar_angle_chest ", planar_angle_chest)
588 # print("Arduino ", highest_axis_index, " ", highest_list_values)
589
590 key = cv2.waitKey(1)
591 if key == ord('h') and help_switch == 0:
592     help_switch = 1
593 elif key == ord('h') and help_switch == 1:
594     help_switch = 0
595
596 # for the user guide
597 cv2.putText(images, 'press h for help', (490, 350), font, 0.5, (0, 0, 255), 1, cv2.LINE_4)
598
599 # draw circle w/ line across
600 images = cv2.circle(images, (600, 75), 36, (0, 0, 0), -1)
601 images = cv2.line(images, (630,55), (578,100), (255,255,255), 2)
602 cv2.putText(images, str(rep_value), (603, 95), font, 0.6, (255, 255, 255), 1, cv2.LINE_4)
603 cv2.putText(images, str(counter), (588, 68), font, 0.6, (255, 255, 255), 1, cv2.LINE_4)
604
605 # for the help screen
606 if help_switch == 1:
607     images = cv2.rectangle(images, (0,0), (640,360), (255, 255, 255), -1)
608
609 cv2.putText(images, 'User guide:', (5,15), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
610 images = cv2.line(images, (5,17), (85,17), (0,0,0), 1)
611 cv2.putText(images, 'Blue',(5,30), font, 0.4, (255, 0, 0), 1, cv2.LINE_4)
612 cv2.putText(images, 'circles represent your', (40,30), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
613 cv2.putText(images, 'skeleton.', (205,30), font, 0.4, (255, 0, 0), 1, cv2.LINE_4)
614 cv2.putText(images, 'The', (5,45), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
615 cv2.putText(images, 'red', (35,45), font, 0.4, (0, 0, 255), 1, cv2.LINE_4)
616 cv2.putText(images, 'lines represent the', (65,45), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
617 cv2.putText(images, 'correct', (205,45), font, 0.4, (0, 0, 255), 1, cv2.LINE_4)
618 cv2.putText(images, 'pose.', (265,45), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
619 cv2.putText(images, 'If the', (5,60), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
620 cv2.putText(images, 'red', (50,60), font, 0.4, (0, 0, 255), 1, cv2.LINE_4)
621 cv2.putText(images, 'lines turn', (80,60), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
622 cv2.putText(images, 'green', (160,60), font, 0.4, (34, 139, 34), 1, cv2.LINE_4)
623 cv2.putText(images, 'correct pose is', (200,60), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
624 cv2.putText(images, 'achieved!', (320,60), font, 0.4, (0, 69, 255), 1, cv2.LINE_4)
625 cv2.putText(images, 'If not, warning messages will appear to tell you how to correct it:', (5, 75), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
626 cv2.putText(images, '- Top left: rotation of arms', (5, 95), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
627 cv2.putText(images, '- Top middle: up or down movement of the arms', (5, 110), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
628 cv2.putText(images, '- Top right: front or back movement of the arms', (5, 125), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
629 cv2.putText(images, '- Bottom left: how far you are from the camera', (5, 140), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
630 cv2.putText(images, '- Bottom middle: elbow corrections', (5, 155), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
631 cv2.putText(images, 'The number at the top right indicates how many times the exercise is completed.', (5, 175), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
632 cv2.putText(images, 'How to:', (5,200), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
633 images = cv2.line(images, (5,202), (55,202), (0,0,0), 1)
634 cv2.putText(images, '1) Go through the gif in the selection screen before proceeding with the exercise.', (5, 215), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
635 cv2.putText(images, '2) Stand 2 to 3m away from the camera.', (5, 230), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
636 cv2.putText(images, '3) The correct pose is shown by the red skeleton in the middle of the screen.', (5, 245), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
637 cv2.putText(images, '4) There is no need to align your skeleton to the red skeleton.', (5, 260), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
638 cv2.putText(images, '5) Follow the red skeleton's pose.', (5, 275), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
639 cv2.putText(images, '6) Any difference between your skeleton and the red skeleton will be shown as', (5, 290), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
640 cv2.putText(images, 'warnings around the screen.', (5, 305), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
641 cv2.putText(images, '7) When one cycle is completed, the counter on the right will increment by one.', (5, 320), font, 0.4, (0, 0, 0), 1, cv2.LINE_4)
642 cv2.putText(images, 'press h to exit', (525, 355), font, 0.4, (0, 0, 255), 1, cv2.LINE_4)
643
644 # make screen full screen
645 cv2.namedWindow("RealSense", cv2.WND_PROP_FULLSCREEN)
646 cv2.setWindowProperty("RealSense",cv2.WND_PROP_FULLSCREEN,cv2.WINDOW_FULLSCREEN)

```

```
648     # Show images
649     cv2.imshow("RealSense", images)
650     # key = cv2.waitKey(1)
651     # Press esc or 'q' to close the image window
652     if key & 0xFF == ord('q') or key == 27:
653         # sys.exit()
654         # cv2.destroyAllWindows()
655         cv2.destroyWindow("RealSense")
656         break
657     elif key == ord('h') and help_switch == 0:
658         help_switch = 1
659     elif key == ord('h') and help_switch == 1:
660         help_switch = 0
661
662
663 finally:
664
665     # Stop streaming
666     pipeline.stop()
667
```

Appendix III (Arduino code)

```
#include <Wire.h>

long accelX, accelY, accelZ;
float gForceX, gForceY, gForceZ;

long gyroX, gyroY, gyroZ;
float rotX, rotY, rotZ;

unsigned long previousMillis = 0;
const long interval = 125;

void setup() {
  Serial.begin(9600);
  Wire.begin();
  setupMPU();
}

void loop() {
  recordAccelRegisters();
  recordGyroRegisters();
  printData();
  delay(100);
}

void setupMPU(){
  Wire.beginTransmission(0b1101000); //This is the I2C address of the MPU (b1101000/b1101001 for AC0 low/high datasheet sec. 9.2)
  Wire.write(0x6B); //Accessing the register 6B - Power Management (Sec. 4.28)
  Wire.write(0b00000000); //Setting SLEEP register to 0. (Required; see Note on p. 9)
  Wire.endTransmission();

  Wire.beginTransmission(0b1101000); //I2C address of the MPU
  Wire.write(0x1B); //Accessing the register 1B - Gyroscope Configuration (Sec. 4.4)
  Wire.write(0x00000000); //Setting the gyro to full scale +/- 250deg./s
  Wire.endTransmission();

  Wire.beginTransmission(0b1101000); //I2C address of the MPU
  Wire.write(0x1C); //Accessing the register 1C - Accelerometer Configuration (Sec. 4.5)
  Wire.write(0b00000000); //Setting the accel to +/- 2g
  Wire.endTransmission();
}

void recordAccelRegisters() {
  Wire.beginTransmission(0b1101000); //I2C address of the MPU
  Wire.write(0x3B); //Starting register for Accel Readings
  Wire.endTransmission();
  Wire.requestFrom(0b1101000,6); //Request Accel Registers (3B - 40)
  while(Wire.available() < 6);
  accelX = Wire.read()<<8|Wire.read(); //Store first two bytes into accelX
  accelY = Wire.read()<<8|Wire.read(); //Store middle two bytes into accelY
  accelZ = Wire.read()<<8|Wire.read(); //Store last two bytes into accelZ
  processAccelData();
}

void processAccelData(){
  gForceX = accelX / 16384.0;
  gForceY = accelY / 16384.0;
  gForceZ = accelZ / 16384.0;
}

void recordGyroRegisters() {
  Wire.beginTransmission(0b1101000); //I2C address of the MPU
```

```
Wire.write(0x43); //Starting register for Gyro Readings
Wire.endTransmission();
Wire.requestFrom(0b1101000,6); //Request Gyro Registers (43 - 48)
while(Wire.available() < 6);
gyroX = Wire.read()<<8|Wire.read(); //Store first two bytes into accelX
gyroY = Wire.read()<<8|Wire.read(); //Store middle two bytes into accelY
gyroZ = Wire.read()<<8|Wire.read(); //Store last two bytes into accelZ
processGyroData();
}

void processGyroData() {
    rotX = gyroX / (131.0*250);
    rotY = gyroY / (131.0*250);
    rotZ = gyroZ / (131.0*250);
}

void printData() {
    recordAccelRegisters();
    recordGyroRegisters();
    Serial.print(gForceX,3);
    Serial.print(",");
    Serial.print(gForceY,3);
    Serial.print(",");
    Serial.print(gForceZ,3);
    Serial.print("\n");
    delay(10);
}
```